

Product Datasheet

AP3M2 Antibody - Azide and BSA Free H00010947-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00010947-B01P

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00010947-B01P



H00010947-B01P

AP3M2 Antibody - Azide and BSA Free

Product Information

Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein G purified
Buffer	PBS (pH 7.4)

Product Description

Host	Mouse
Gene ID	10947
Gene Symbol	AP3M2
Species	Human
Specificity/Sensitivity	This product is specific for Human AP3M2 purified MaxPab mouse polyclonal antibody (B01P) [Gene ID: 10947].
Immunogen	AP3M2 (NP_006794.1, 1 a.a. - 418 a.a.) full-length human protein. MIHSLFLINSSGDIFLEKHWKSVVSRVCDYFFEAQERATEAENVPPVIPTPHHY LLSVYRHKIFFVAVIQTEVPPLFVIEFLHRVVDTFQDYFGVCSEPVIKDNVVVVYE VLEEMLDNGFPLATESNILKELIKPPTILRTVVNTITGSTNVGDQLPTGQLSVVP WRRTGVKYTNNEAYFDVIEEIDAIDKSGSTITAEIQGVIDACVKLTGMPDLTSLF MNPRLDDVSFHPCVRFKRWESERILSFIPPDGNFRLLSYHVS AQNLVAIPVYV KHNISFRDSSSLGRFEITVGPKQTMGKTIEGVTVT SQMPKGV LNM SLTPSQGT HTFDPVTKMLSVDVGKINPQKLPSLKG TMSLQAGASKPDENPTINLQFKIQQLA ISGLKVNRLDMYGEKYKPFKG I KYMTKAGKFQVRT
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Mouse polyclonal antibody raised against a full-length human AP3M2 protein.

Images

Western Blot: AP3M2 Antibody [H00010947-B01P] - Analysis of AP3M2 expression in transfected 293T cell line by AP3M2 polyclonal antibody. Lane 1: AP3M2 transfected lysate(45.98 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to H00010947-B01P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
H00010947-P01-10ug	Recombinant Human AP3M2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00010947-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

