

Product Datasheet

MIG2/Kindlin-2 Antibody (2G11) - Azide and BSA Free H00010979-M09

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00010979-M09

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00010979-M09



H00010979-M09

MIG2/Kindlin-2 Antibody (2G11) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2G11
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4

Product Description	
Description	Novus Biologicals Mouse MIG2/Kindlin-2 Antibody (2G11) - Azide and BSA Free (H00010979-M09) is a monoclonal antibody validated for use in WB, ELISA and ICC/IF. Anti-MIG2/Kindlin-2 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	10979
Gene Symbol	FERMT2
Species	Human, Mouse
Reactivity Notes	Human. Other species not tested.
Specificity/Sensitivity	PLEKHC1 (2G11)
Immunogen	PLEKHC1 (AAH17327, 1 a.a. ~ 680 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MALDGI RMPDGCYADGTWELSVHVTDLNRDVTLRVTGEVHIGGVMLKLVEKLD VKKDWS DHALWWEKKRTWLLKTHWTLDKYGIQADAKLQFTPQHKLLRLQLPN MKYVKVKVNFSDRVFKA VSDICKTFNIRHPEELSLLKKPRDPTKKKKKKLDDQS EDEALELEGPLITPGSGSIYSSPGLYSKTMPTTYDAHDGSPLSPTS AWFGDSAL SEGNPGILAVSQPITSPEILAKMFKPQALLDKAKINQGWLDSSRSLMEQDVKEN EALLRKFYYSFFDLNPKYDAIRINQLYEQAKWAILLEEIECTEEEMMMFAALQY HINKLSIMTSENHLNNSDKEVDEVDAALSDLEITLEGGKTSTILGDITSIPELADYI KVFKPKKLT LKGYKQYWCTFKDTSISCYKSKEESSGTPAHQMNLRGCEVTPDV NISGQKFNIKLLIPVAEGMNEIWLRC DNEKQYAHWMAACRLASKGKTMADSSY NLEVQNILSS LKMQHLPDPQLIPEQITTDITPECLVSPRYLKKYKNKQITARILEA HQNV AQMSLIEAKMRFIQAWQSLPEFGITHFIARFQGGKKEELIGIAYNRLIRMD ASTGDAIKTWRF SNMKQWNVNWEIKMVTVEFADEVRLSFICTEVDCKVVHEFI GGYIFLSTRAKDQNESLDEEMFYKLTSGWV
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Sandwich ELISA
Recommended Dilutions	Western Blot 1:500, ELISA, Immunocytochemistry/ Immunofluorescence, Sandwich ELISA




Application Notes


Antibody reactive against recombinant protein on ELISA. It has also been used for Western Blot.

Images


Western Blot: MIG2/Kindlin-2 Antibody (2G11) [H00010979-M09] - Analysis of PLEKHC1 expression in transfected 293T cell line by PLEKHC1 monoclonal antibody (M09), clone 2G11. Lane 1: PLEKHC1 transfected lysate (77.9 KDa). Lane 2: Non-transfected lysate.



Immunocytochemistry/Immunofluorescence: MIG2/Kindlin-2 Antibody (2G11) [H00010979-M09] - Analysis of monoclonal antibody to FERMT2 on HeLa cell. Antibody concentration 10 ug/ml



Sandwich ELISA: MIG2/Kindlin-2 Antibody (2G11) [H00010979-M09] - Detection limit for recombinant GST tagged FERMT2 is 1 ng/ml as a capture antibody.

**Publications**

Halbert D, Domenyuk V, Spetzler D et al. Aptamers and uses thereof United States Patent Application US 9958448 B2 2018-01-01

Dowling JJ, Vreede AP, Kim S et al. Kindlin-2 is required for myocyte elongation and is essential for myogenesis. BMC Cell Biol. 2008-07-08 [PMID: 18611274]

**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00010979-M09

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
H00010979-P01-10ug	Recombinant Human MIG2/Kindlin-2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00010979-M09

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

