

Product Datasheet

SF3B2 Antibody (5D2) - Azide and BSA Free H00010992-M01

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00010992-M01

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00010992-M01



H00010992-M01

SF3B2 Antibody (5D2) - Azide and BSA Free

Product Information

Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	5D2
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4

Product Description

Host	Mouse
Gene ID	10992
Gene Symbol	SF3B2
Species	Human, Mouse, Rat
Specificity/Sensitivity	SF3B2 (5D2)
Immunogen	SF3B2 (NP_006833, 592 a.a. ~ 645 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. YEGKEFETRLKEKKPGDLSDELRLSLGMPVGPNAHKVPPPWLIAMQRYGPPPS Y
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details


Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin
Application Notes	Antibody reactive against cell lysate and recombinant protein for Western Blot. Has also been used for immunofluorescence, immunohistochemistry (paraffin), and ELISA.

Images


Western Blot: SF3B2 Antibody (5D2) [H00010992-M01] - SF3B2 monoclonal antibody (M01), clone 5D2 Analysis of SF3B2 expression in Hela S3 NE.




Immunocytochemistry/Immunofluorescence: SF3B2 Antibody (5D2) [H00010992-M01] - Analysis of monoclonal antibody to SF3B2 on HeLa cell. Antibody concentration 60 ug/ml.




Immunohistochemistry-Paraffin: SF3B2 Antibody (5D2) [H00010992-M01] - Analysis of monoclonal antibody to SF3B2 on formalin-fixed paraffin-embedded human kidney. Antibody concentration 6 ug/ml.



Western Blot: SF3B2 Antibody (5D2) [H00010992-M01] - SF3B2 monoclonal antibody (M01), clone 5D2. Analysis of SF3B2 expression in PC-12.



Western Blot: SF3B2 Antibody (5D2) [H00010992-M01] - SF3B2 monoclonal antibody (M01), clone 5D2. Analysis of SF3B2 expression in NIH/3T3.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00010992-M01

NB800-PC9	HeLa Nuclear Cell Lysate
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00010992-M01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

