

Product Datasheet

KIF3A Antibody - Azide and BSA Free H00011127-B01P-50ug

Unit Size: 50 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00011127-B01P-50ug

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00011127-B01P-50ug



H00011127-B01P-50ug

KIF3A Antibody - Azide and BSA Free

Product Information	
Unit Size	50 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Novus Biologicals Mouse KIF3A Antibody - Azide and BSA Free (H00011127-B01P-50ug) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	11127
Gene Symbol	KIF3A
Species	Human, Rat
Immunogen	KIF3A (AAH45542.1, 1 a.a. - 702 a.a.) full-length human protein. MPINKSEKPESCDNVKKVVRCRPLNEREKSMCYKQAVSVDEMGRGTITVHKTD SNEPPKTFTFDVFPGESKQLDVYNLTARPIIDSVLEGYNGTIFAYGQTGTGKTF TMEGVRAIPELRGIIPNSFAHIFGHIKAEGDTRFLVRVSYLGIYNEEVRDLLGKD QTQRLEVIERPDVGVIKDL SAYVVNNADDMDRIMTLGHKNRSGATNMNEHS SRSHAIFTITIECSEKIDGNMHVRMGKLHLVDLAGSERQAKTGATGQRLKEAT KINLSLSTLGNVISALVDGKSTHVPYRNSKLTRLLQDSLGGNSKTMMCANIGPA DNYDETISTLRYANRAKNIKNKARINEDPKDALLRQFQKEIEELKKKLEEGEEIS GSDISGSEEDDDEEGEVGEDGEKRKKRRDQTGKKKVSPDKMIEMQAKIDEER KALETKLDMEEEERNKARAELEKREKDLLKAQQEHQSLEKLSALEKKVIVGGV DLLAKAEEQEKLLEESNMELEERRKRAEQLRRELEEKEQERLDIEEKYTSLQEE AQGKTKKLKKVWTMLMAAKSEMADLQQEHQREIEGLLENIRQLSRELRLQMLII DNFIPRDYQEMIENYVHWNE DIGEWQLKCVAYTGNNMRKQTPVPDKKEKDPF EVDLSHVYLAYTEESLRQSLMKLERPRTSKGKARPKTGRRKRSAPETVIDSL Q
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot




Images

Western Blot: KIF3A Antibody [H00011127-B01P-50ug] - MaxPab polyclonal antibody. Western Blot analysis of KIF3A expression in rat brain.



Western Blot: KIF3A Antibody [H00011127-B01P-50ug] - Analysis of KIF3A expression in transfected 293T cell line (H00011127-T01) by KIF3A MaxPab polyclonal antibody. Lane1:KIF3A transfected lysate (77.22 KDa). Lane2:Non-transfected lysate.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00011127-B01P-50ug

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
H00011127-P01-10ug	Recombinant Human KIF3A GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00011127-B01P-50ug

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

