

Product Datasheet

IRAK3 Antibody (1C9) - Azide and BSA Free H00011213-M07

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00011213-M07

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00011213-M07



H00011213-M07


IRAK3 Antibody (1C9) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1C9
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	Protein A purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Host	Mouse
Gene ID	11213
Gene Symbol	IRAK3
Species	Human, Mouse, Rat
Specificity/Sensitivity	Reacts with interleukin-1 receptor-associated kinase 3.
Immunogen	IRAK3 (AAH57800, 1 a.a. ~ 596 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MAGNWGARGALSAHTLLFDLPPALLGELCAVLDS CDGALGWRGLAERLSSSR LDVRHIEKYVDQ GKSGTRELLWSWAQKNKTIGDLLQVLQEMGHRRAIHLITNY GAVLSPSEKSYQEGGFPNILFKETANVTVDNVLIP EHNKGVLLKSSISFQNIIEG TRNFHKDFLIGEGEIFEVYRVEIQNLTYAVKLFKQEKKMQCKKHWKRFLSELEA LLLFHHPSILELAAYFTETEFCL IYPYMRNGTLFDRLQCVGDTAPLPWHIRIGILI GISKAIHYLHNVQPCSVICGSISSANILLDDQFQPKLTD FAMAHRSHLEHQSCTI NMTSSSSKHLWYMPEEYIRQGKLSIKTDVYSFGIVIMEVLTGCRVVLDDPKHIQL RDLLRELMEKRGLDSCLSFLDKKVP PCPRNFSAKLFCLAGRCAATRAKLRPSM DEVLNTLESTQASLYFAEDPPTSLKSFRCPSPLFLENVPSIPVEDDESQNNLL PSDEGLRIDRMTQKTPFECSQSEVMFLSLDKKPESKRNEEACNMPSSSCEES WFPKYIVPSQDLRPYKVNIDPSSEAPGHSCRSRPVESSCSSKFSWDEYEYKK E
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot, ELISA
Application Notes	This antibody is reactive against cell lysate in western blot and recombinant protein in ELISA.




Images


Western Blot: IRAK3 Antibody (1C9) [H00011213-M07] - Analysis of IRAK3 expression in Raw 264.7.




Western Blot: IRAK3 Antibody (1C9) [H00011213-M07] - Analysis of IRAK3 expression in HeLa.



Western Blot: IRAK3 Antibody (1C9) [H00011213-M07] - Analysis of IRAK3 expression in PC-12.



Western Blot: IRAK3 Antibody (1C9) [H00011213-M07] - Analysis of IRAK3 expression in Jurkat (Cat # L017V1).



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00011213-M07

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
H00011213-P01-10ug	Recombinant Human IRAK3 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00011213-M07

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

