

Product Datasheet

DIS3 Antibody - Azide and BSA Free H00022894-B01P

Unit Size: 50 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 2 Publications: 4

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00022894-B01P

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00022894-B01P



H00022894-B01P

DIS3 Antibody - Azide and BSA Free

Product Information	
Unit Size	50 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein G purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Novus Biologicals Mouse DIS3 Antibody - Azide and BSA Free (H00022894-B01P) is a polyclonal antibody validated for use in WB and ICC/IF. Anti-DIS3 Antibody: Cited in 4 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	22894
Gene Symbol	DIS3
Species	Human, Mouse, Chicken
Immunogen	KIAA1008 (AAH56143.1, 1 a.a. - 958 a.a.) full-length human protein. MLKSKTFLKKTRAGGVMKIVREHYLRDDIGCGAPGCAACGGAHEGPALEPQP QDPASSVCPQPHYLLPDTNVLLHQIDVLEDPAIRNVIVLQTVLQEVNRNSAPVY KRIRDVTNNQEKHFYFTTNEHHRETYVEQEQQENANDRNDRAIRVAAKWYNE HLKKMSADNQLQVIFITNDRRNKEKAIEEGIPAFTCEEYVKSLTANPELIDRLACL SEEGNEIESGKIIFSEHLPLSKLQQGIKSGTYLQGTFRASRENYLEATVWIHGDN EENKEIILQGLKHLNRAVHEDIVAVELLPKSQWVAPSSSVVLHDEGQNEEDVEKE EERERMLKTAVSEKMLKPTGRVVGIIKRNWRPYCGMLSLSKSDIKESRRHLFTPA DKRIPRIETRQASTLEGRRIIVADGWPRNSRYPNGHFVRNLGVDVGEKETETE VLLLEHDVPHQPFSQAVLSFLPKMPWSITEKDMKNREDLRHLCICSVDPPGCT DIDDALHCRELENGNLEVGVIADVSHFIRPGNALDQESARRGTTVYLCEKRID MYPELLSSNLCSLKCDVDRLAFSCIWEMNHNAEILKTKFTKSVINSKASLTIAEA QLRIDSANMDDITSLRGLNKLAKILKKRRIEKGALTSSPEVRFHMDSETHDPI DLQTKELRETNMVEEFMLLANISVAKKIEEFSEHALLRKHPAPPPSNYEILVK AARSRNLEIKTDTAKSLAESLDQAESPTFPYLNTLLRILATRCMMQAVYFCSGM DNDFHHYGLASPIYTHFTSPIRRYADVIVHRLLAIVAIGADCTYPELTDKHKLADIC KNLNFRHKMAQYAQRASVAFHTQLFFKSKGIVSEEAYILFVRKNAIVVLIPKYGL EGTVFFEEKDKPNPQLIYDDEIPSLKIEDTVFHVFDKVKVKIMLDSSNLQHQKIR MSLVEPQIPGISIPTDTSNMDLNGPKKKKMKLGK
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:100-1:2000, Immunocytochemistry/ Immunofluorescence 1:10-1:500
Application Notes	This antibody is useful for Western Blot, Immunofluorescence




Images

Western Blot: DIS3 Antibody [H00022894-B01P] - analysis of DIS3 in HeLa whole cell extracts treated with (1) si-Control and (2) DIS3 siRNA using anti-DIS3 antibody. Image from verified customer review.



Western Blot: DIS3 Antibody [H00022894-B01P] - Analysis of DIS3 expression in transfected 293T cell line by DIS3 polyclonal antibody. Lane 1: KIAA1008 transfected lysate(105.38 KDa). Lane 2: Non-transfected lysate.



Immunocytochemistry/Immunofluorescence: DIS3 Antibody [H00022894-B01P] - Analysis of purified antibody to KIAA1008 on HeLa cell. (antibody concentration 10 ug/ml)



Publications

Scott, Micha M, Arpan R et al. Feedback from nuclear RNA on transcription promotes robust RNA concentration homeostasis in human cells. *Cell Syst.* 2022-05-18 [PMID: 35613616]

Ogami K, Richard P, Chen Y et al. An Mtr4/ZFC3H1 complex facilitates turnover of unstable nuclear RNAs to prevent their cytoplasmic transport and global translational repression. *Gens & Dev.* 2017-07-21 [PMID: 28733371]

Xie M, Li M, Vilborg A et al. Mammalian 5'-Capped MicroRNA Precursors that Generate a Single MicroRNA. *Cell.* 2013-12-19 [PMID: 24360278]

Hsin JP, Li W, Hoque M et al. RNAP II CTD tyrosine 1 performs diverse functions in vertebrate cells. *Elife* (Cambridge) 2014-06-04 [PMID: 24842995] (WB, Chicken)

Details:

EXOSC3 antibody used for WB on lysates of 26r (DT40 cells derived Rpb1 derivative containing a CTD with 26 YSPTSPS repeats) and 25F+Y cells (Rpb1-Y1F derivative in which only a single F, in the C terminal-most heptad, was changed back to Y). WB data shown in Supplement Figure 4.

**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00022894-B01P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
H00022894-Q01-10ug	Recombinant Human DIS3 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00022894-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

