

Product Datasheet

Senataxin Antibody (2A9) - Azide and BSA Free H00023064-M05

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00023064-M05

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00023064-M05



H00023064-M05

Senataxin Antibody (2A9) - Azide and BSA Free

Product Information

Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2A9
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	Protein A purified
Buffer	In 1x PBS, pH 7.4

Product Description


Description	Novus Biologicals Mouse Senataxin Antibody (2A9) - Azide and BSA Free (H00023064-M05) is a monoclonal antibody validated for use in WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	23064
Gene Symbol	SETX
Species	Human, Mouse, Rat
Specificity/Sensitivity	SETX (2A9)
Immunogen	SETX (NP_055861.2, 2579 a.a. ~ 2676 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. VVHQDLSHIQQPAAVVAALSSHKPPVRGEPPAASPEASTCQSKCDDPEEELCH RREARAFSEGEQEKCGETHHTRNSRWDRKRTLEQEDSSSKRKL
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details


Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:500, ELISA, Immunocytochemistry/ Immunofluorescence
Application Notes	Antibody reactive against cell lysate and recombinant protein for Western Blot. Has also been used for immunofluorescence and ELISA.

Images


Western Blot: Senataxin Antibody (2A9) [H00023064-M05] - Analysis of SETX expression in PC-12 (Cat # L012V1).




Immunocytochemistry/Immunofluorescence: Senataxin Antibody (2A9) [H00023064-M05] - Analysis of monoclonal antibody to SETX on HeLa cell. Antibody concentration 10 ug/ml




Western Blot: Senataxin Antibody (2A9) [H00023064-M05] - SETX expression in Raw 264.7



Western Blot: Senataxin Antibody (2A9) [H00023064-M05] - Analysis of SETX expression in NIH/3T3 (Cat # L018V1).



ELISA: Senataxin Antibody (2A9) [H00023064-M05] - Detection limit for recombinant GST tagged SETX is 0.3 ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00023064-M05

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00023064-M05

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

