

Product Datasheet

CDC42EP4 Antibody (3G10) - Azide and BSA Free H00023580-M05

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00023580-M05

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00023580-M05



H00023580-M05

CDC42EP4 Antibody (3G10) - Azide and BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3G10
Preservative	No Preservative
Isotype	IgG2b Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4

Product Description	
Description	Novus Biologicals Mouse CDC42EP4 Antibody (3G10) - Azide and BSA Free (H00023580-M05) is a monoclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	23580
Gene Symbol	CDC42EP4
Species	Human
Specificity/Sensitivity	CDC42EP4 - CDC42 effector protein (Rho GTPase binding) 4 (3G10)
Immunogen	CDC42EP4 (NP_036253, 163 a.a. ~ 225 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. VPRRNGAAGPHSPDLLDEQAFGLTDLPPVVPKATYGLKHAESIMSFHIDLGPS MLGDVLSIM
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:500, ELISA
Application Notes	Antibody Reactive Against Recombinant Protein with GST tag on ELISA and Western Blot. GST tag alone is used as a negative control.

Images

Western Blot: CDC42EP4 Antibody (3G10) [H00023580-M05] - Analysis of CDC42EP4 expression in transfected 293T cell line by CDC42EP4 monoclonal antibody (M05), clone 3G10. Lane 1: CDC42EP4 transfected lysate (Predicted MW: 38 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to H00023580-M05

NBP2-33376H

Blue Marker Antibody (6F4-F6) [HRP]

HAF007

Goat anti-Mouse IgG Secondary Antibody [HRP]

NB7539

Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]

NBP1-43317-0.5mg

Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00023580-M05

Earn gift cards/discounts by submitting a publication using this product:

www.novusbio.com/publications

