

Product Datasheet

IGSF4A/SynCAM1/CADM1 Antibody - Azide and BSA Free H00023705-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00023705-D01P

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00023705-D01P



H00023705-D01P

IGSF4A/SynCAM1/CADM1 Antibody - Azide and BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)

Product Description	
Host	Rabbit
Gene ID	23705
Gene Symbol	CADM1
Species	Human
Specificity/Sensitivity	Reacts with cell adhesion molecule 1.
Immunogen	CADM1 (AAH47021.1, 1 a.a. - 387 a.a.) full-length human protein. MASVVLPSGSQCAAAAAAAAAAPPGLRLRLLLLLFSAAALIPTGDGQNLFTKDVTVI EGEVATISCQVNKSDDSVIQLLNPNRQTIYFRDFRPLKDSRFQLLNFSSELKVS LTNVSISDEGRYFCQLYTDPPQESYTTITVLVPPRNLMIDIQKDTAVEGEEIEVN CTAMASKPATTIRWFKGNTELKKGKSEVEEWSDMYTVTSQLMLKVHKEDDGVP VICQVEHPAVTGNLQTQRYLEVQYKPQVHIQMTYPLQGLTREGDALELTCEAIG KPQPMVMTWVRVDDEMPQHAVLSGPNLFINNLNKTNGTYRCEASNIVGKAH SDYMLYVYDSRAGEEGSIRAVDHA VIGGVVAVVVVFAMLCLLILGRYFARHKGL FSLTSSPRIK
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	This antibody is reactive against transfected lysate in western blot, and as a detection antibody in ELISA.

Images

Western Blot: IGSF4A/SynCAM1/CADM1 Antibody [H00023705-D01P] - Analysis of CADM1 expression in transfected 293T cell line by CADM1 polyclonal antibody. Lane 1: CADM1 transfected lysate(42.70 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to H00023705-D01P

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NB300-186PEP	IGSF4A/SynCAM1/CADM1 Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00023705-D01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

