

Product Datasheet

PHGDH Antibody (4A3-1D6) - Azide and BSA Free H00026227-M01

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00026227-M01

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00026227-M01



H00026227-M01

PHGDH Antibody (4A3-1D6) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	4A3-1D6
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Novus Biologicals Mouse PHGDH Antibody (4A3-1D6) - Azide and BSA Free (H00026227-M01) is a monoclonal antibody validated for use in IHC, WB, ELISA, ICC/IF and IP. Anti-PHGDH Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	26227
Gene Symbol	PHGDH
Species	Human
Reactivity Notes	Human. Other species not tested.
Specificity/Sensitivity	PHGDH - phosphoglycerate dehydrogenase
Immunogen	PHGDH (AAH11262.1, 1 a.a. ~ 533 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MAFANLRKVLISDSLDPCCRKILQDGGQLQVVEKQNLKEELIAELQDCEGLIVRS ATKVTADVINA AEKLQVVG RAGTGV DNV DLEAATRKGILVMNTPNGNSLSAAEL TCGMIMCLARQIPQATASMKGDKWERKKFMGTENGLTLGILGLGRIGREVAT RMQSFGMKTIGYDPIISPEVSASFGVQQLPLEEIWPLCDFITVHTPLLPSTTGLL NDNTFAQCKKGVRV VNCARGGIVDEGALLRALQSGQCAGAALDVFTTEPPRD RALVDHENVISCPHLGASTKEAQSRCGEEIAVQFVDMVKGKSLTGVVNAQALT SAFSPHTKPWIGLAEALGTL MRAWAGSPKGTIQVITQGTSLKNAGNCLSPAVIV GLLKEASKQADVNLVNAKLLVKEAGLNVTTSHSPAAPGEQGFGECLLAVALAG APYQAVGLVQGTT PVLQGLNGAVFRPEVPLRRDLPLLLFRTQTSDPAMLPTMI GLLAEAGVRLLSYQTSLVSDGETWHVMGISSLLPSLEAWKQHVTEAFQFHF
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Immunohistochemistry-Paraffin
Application Notes	Antibody reactive against recombinant protein, tissue lysate and cell lysate for Western Blot. Has also been used for immunofluorescence, immunohistochemistry (paraffin), and ELISA.




Images


Western Blot: PHGDH Antibody (4A3-1D6) [H00026227-M01] - PHGDH monoclonal antibody (M01), clone 4A3-1D6. Analysis of PHGDH expression in human liver.




Immunocytochemistry/Immunofluorescence: PHGDH Antibody (4A3-1D6) [H00026227-M01] - Analysis of monoclonal antibody to PHGDH on HeLa cell. Antibody concentration 1 ~ 10 ug/ml.




Immunohistochemistry-Paraffin: PHGDH Antibody (4A3-1D6) [H00026227-M01] - Analysis of monoclonal antibody to PHGDH on formalin-fixed paraffin-embedded human prostate. Antibody concentration 3 ug/ml.




Western Blot: PHGDH Antibody (4A3-1D6) [H00026227-M01] - PHGDH monoclonal antibody (M01), clone 4A3-1D6. Analysis of PHGDH expression in A-431.




Western Blot: PHGDH Antibody (4A3-1D6) [H00026227-M01] - PHGDH monoclonal antibody (M01), clone 4A3-1D6 Analysis of PHGDH expression in Jurkat.




Western Blot: PHGDH Antibody (4A3-1D6) [H00026227-M01] - Analysis of PHGDH expression in transfected 293T cell line by PHGDH monoclonal antibody (M01), clone 4A3-1D6. Lane 1: PHGDH transfected lysate (56.7 kDa). Lane 2: Non-transfected lysate.



Immunohistochemistry-Paraffin: PHGDH Antibody (4A3-1D6) [H00026227-M01] - Analysis of monoclonal antibody to PHGDH on formalin-fixed paraffin-embedded human lymph node. Antibody concentration 1 ~ 10 ug/ml.



Immunoprecipitation: PHGDH Antibody (4A3-1D6) [H00026227-M01] - Analysis of PHGDH transfected lysate using anti-PHGDH monoclonal antibody and Protein A Magnetic Bead, and immunoblotted with PHGDH MaxPab rabbit polyclonal antibody.



Publications

Halbert D, Domenyuk V, Spetzler D et al. Aptamers and uses thereof United States Patent Application US 9958448 B2 2018-01-01

Pollari S, Kakonen SM, Edgren H et al. Enhanced serine production by bone metastatic breast cancer cells stimulates osteoclastogenesis. Breast Cancer Res Treat. 2010-03-30 [PMID: 20352489]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to H00026227-M01

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
NBP1-87311PEP	PHGDH Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00026227-M01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

