

Product Datasheet

BRPF3 Antibody - Azide and BSA Free H00027154-B01P

Unit Size: 50 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00027154-B01P

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00027154-B01P



H00027154-B01P

BRPF3 Antibody - Azide and BSA Free

| Product Information | |
|---------------------|--|
| Unit Size | 50 ug |
| Concentration | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | No Preservative |
| Isotype | IgG |
| Purity | Protein G purified |
| Buffer | PBS (pH 7.4) |

| Product Description | |
|---------------------|---|
| Description | Novus Biologicals Mouse BRPF3 Antibody - Azide and BSA Free (H00027154-B01P) is a polyclonal antibody validated for use in WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Mouse |
| Gene ID | 27154 |
| Gene Symbol | BRPF3 |
| Species | Human |
| Reactivity Notes | This product is reactive against Human. |
| Immunogen | BRPF3 (NP_056510.1, 1 a.a. - 1205 a.a.) full-length human protein. MRKPRRKSQRNAEGRSPSPYSLKCSPTRETLTYAQAQRIVEVDIDGRLHRISIDPLKIITEDELTAQDITECNSNKENSEQPQFPGKSKKPKSSKGGKKKESCSKHASGTSFHLPPQPSFRMVDSGIQPEAPPLPAAYYRYIEKPPEDLDAEVEYDMDEEDLAWLDMVNEKRRVDGHSLSVSADTFELLVDRLEKESYLESRSSGAQQSLIDEDAFCCVCLDDECHNSNVILFCDICNLAVHQECYGVPIPEGQWLCRCCLQSPSRPVDICILCPNKGGAFFKQTSBGHWAHVVCIAWIPEVCFANTVFLEPIEGIDNIPPARWKLTCYICKQKGLGAAIQCHKVNCYTAFHVTCARAGLFMKIEPMRETSLNGTIFTVRKTAYCEAHSPPGAATARRKGDSPRSISETGDEEGLKEGDGEEEEEEEEVEE EEQEAQGGVSGSLKGVPKKSKMSLQKQIKKEPEEAGQDTPSTLPLAVPQIPS YRLNKICSGLSFQRKNQFMQRLHNYWLLKRQARNGVPLIRRLHSHLQSQRNAE QREQDEKTSVKEELKYWQKLRHDLERARLLIELIRKREKLKREQVKVQQAAM EELMPFNVLLRRTLLDQLQKDPAHIFAEPVNLSEVPDYLEFISKPMDFSTMRRL LESHLYRTLEEFEDFNLIVTNCMKYNAKDTIFHRAAVRLRDLGGAILRHARRQA ENIGYDPERGTHLPESPKLEDFYRFSWEDVDNIPENRAHLSPEVQLKELLEKL DLVSAMRSSGARTRRVLLRREINALRQKLAQPPPPQPPSLNKTVSNGELPAG PQGDAVLEQALQEEPEDDGDGDRDDSKLPPPPTLEPTGPAPSLSEQESPEPPT LKPINDSKPPSRFLKPRKVEEDELLEKSPQLQNGNEPLQRLSDNGINRLSLMAP DTPAGTPLSGVGRRTSVLFKAKNGVKLQRSPDRVLENGEDHGVAGSPASPA SIEEERHSRKRPRSRSCSESEGERSPQQEEETGMTNGFGKHTESGSDSECSL GLSGGLAFEACSGLTPPKRSRGKPALSRVPFLEGVNGSDYNGSGRSLLLPFE DRGDLEPLELVWAKCRGYPSYPALIIDPKMPREGLLHNGVPIPVPLDVLKLG QKQAEAGEKFLVLFNDKRTWQWLPDRDKVPLGVEDTVDKLKMLEGRKTSIR KSVQVAYDRAMIHLSRVRGPHSFVTSSYL |
| Notes | This product is produced by and distributed for Abnova, a company based in Taiwan. |


| Product Application Details | |
|-----------------------------|---|
| Applications | Western Blot, Immunocytochemistry/ Immunofluorescence |




| | |
|------------------------------|--|
| Recommended Dilutions | Western Blot, Immunocytochemistry/ Immunofluorescence |
| Application Notes | This antibody is useful for Western Blot, Immunofluorescence |

Images

Western Blot: BRPF3 Antibody [H00027154-B01P] - Analysis of BRPF3 expression in transfected 293T cell line by BRPF3 polyclonal antibody. Lane 1: BRPF3 transfected lysate(132.55 KDa). Lane 2: Non-transfected lysate.



Immunocytochemistry/Immunofluorescence: BRPF3 Antibody [H00027154-B01P] - Analysis of purified antibody to BRPF3 on HeLa cell. (antibody concentration 10 ug/ml)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00027154-B01P

| | |
|----------------|--|
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP] |
| HAF007 | Goat anti-Mouse IgG Secondary Antibody [HRP] |
| NB7539 | Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP] |
| NBP1-97019-5mg | Mouse IgG Isotype Control |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00027154-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

