

Product Datasheet

PD-L1 Antibody Pair H00029126-AP41

Unit Size: 1 Set

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00029126-AP41

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00029126-AP41



H00029126-AP41

PD-L1 Antibody Pair


Product Information	
Unit Size	1 Set
Concentration	Concentration of individual antibodies may be found on the vial label. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Buffer	1x PBS, pH 7.4

Product Description	
Description	Quality control test: Standard curve using recombinant protein (H00029126-Q01) as an analyte. This antibody pair set comes with matched antibody pair to detect and quantify protein level of human CD274.
Gene ID	29126
Gene Symbol	CD274
Species	Human
Kit Components	Capture antibody: mouse monoclonal anti-CD274 IgG2b Kappa (100 ug), Detection antibody: biotinylated mouse monoclonal anti-CD274 IgG2b Kappa (50 ug)
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Standard Curve Range	0.3 ng/ml to 100 ng/ml
Assay Type	Sandwich ELISA

Product Application Details	
Applications	Sandwich ELISA
Recommended Dilutions	Sandwich ELISA 0.3 ng/ml - 100 ng/ml
Application Notes	Reagents are sufficient for at least 3-5 x 96 well plates using recommended protocols.

Images

Sandwich ELISA: CD274 Antibody Pair [H00029126-AP41] - Detection sensitivity ranging from 0.3 ng/ml to 100 ng/ml.

**Publications**

Ingersoll SA, Laval J, Forrest OA et al. Mature cystic fibrosis airway neutrophils suppress T cell function: evidence for a role of arginase 1 but not programmed death-ligand 1. *J Immunol.* 2015-04-29 [PMID: 25926674]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00029126-AP41

NBP1-76769PEP-0.1mg	PD-L1 Antibody Blocking Peptide
210-TA-005	TNF-alpha [Unconjugated]
NBP1-76769	PD-L1 Antibody - BSA Free
6507-IL-010/CF	IL-4 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Antibody Pairs are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00029126-AP41

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

