

# Product Datasheet

## KChIP2 Antibody - Azide and BSA Free H00030819-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00030819-D01P](http://www.novusbio.com/H00030819-D01P)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00030819-D01P](http://www.novusbio.com/reviews/destination/H00030819-D01P)



**H00030819-D01P**


KChIP2 Antibody - Azide and BSA Free

| <b>Product Information</b>         |   |
|------------------------------------|---|
| <b>Unit Size</b>                   | 0.1 mg  |
| <b>Concentration</b>               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.  |
| <b>Storage</b>                     | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.  |
| <b>Clonality</b>                   | Polyclonal  |
| <b>Preservative</b>                | No Preservative   |
| <b>Isotype</b>                     | IgG   |
| <b>Purity</b>                      | Protein A purified  |
| <b>Buffer</b>                      | PBS (pH 7.4)  |
| <b>Product Description</b>         |   |
| <b>Description</b>                 | Novus Biologicals Rabbit KChIP2 Antibody - Azide and BSA Free (H00030819-D01P) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.  |
| <b>Host</b>                        | Rabbit  |
| <b>Gene ID</b>                     | 30819   |
| <b>Gene Symbol</b>                 | KCNIP2  |
| <b>Species</b>                     | Human   |
| <b>Specificity/Sensitivity</b>     | KCNIP2 - Kv channel interacting protein 2,  |
| <b>Immunogen</b>                   | KCNIP2 (NP_775283.1, 1 a.a. - 270 a.a.) full-length human protein.<br>MRGQGRKESLSDSRDLGSDYDQLTGHPGPTKKALKQRFLKLLPCCGPQALP<br>SVSETLAAPASLRPHRPRLLDPDSVDDEFELSTVCHRPEGLEQLQEQTKEFTRK<br>ELQVLYRGGFKNECPGIVNEENFKQIYSQFFPQGDSSTYATFLFNAFDTNHDGS<br>VSFEDFVAGLSVILRGTVDDRLNWFNLYDLNKDGCITKEEMLDIMKSIYDMMG<br>KYTYPALREEAPREHVESFFQKMDRKNKDGVVTTIEEFIESCQKDENIMRSMQLF<br>DNVI |
| <b>Notes</b>                       | This product is produced by and distributed for Abnova, a company based in Taiwan.  |
| <b>Product Application Details</b> |   |
| <b>Applications</b>                | Western Blot  |
| <b>Recommended Dilutions</b>       | Western Blot 1:500  |
| <b>Application Notes</b>           | Antibody reactive against tissue lysate and transfected lysate for western blot. It has also been used for ELISA.   |



**Images**


Western Blot: KChIP2 Antibody [H00030819-D01P] - Analysis of KCNIP2 expression in human kidney.



Western Blot: KChIP2 Antibody [H00030819-D01P] - Analysis of KCNIP2 expression in A-431.



Western Blot: KChIP2 Antibody [H00030819-D01P] - Analysis of KCNIP2 expression in transfected 293T cell line by KCNIP2 polyclonal antibody. Lane 1: KCNIP2 transfected lysate(30.90 KDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00030819-D01P**

---

|             |   |
|-------------|---|
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP]                 |
| HAF008      | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160      | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891  | Rabbit IgG Isotype Control                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00030819-D01P](http://www.novusbio.com/reviews/submit/H00030819-D01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

