

Product Datasheet

RNF111 Antibody - Azide and BSA Free H00054778-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00054778-B01P

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00054778-B01P



H00054778-B01P

RNF111 Antibody - Azide and BSA Free

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)

Product Description	
Description	Novus Biologicals Mouse RNF111 Antibody - Azide and BSA Free (H00054778-B01P) is a polyclonal antibody validated for use in IHC and WB. Anti-RNF111 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	54778
Gene Symbol	RNF111
Species	Human, Mouse, Rat
Specificity/Sensitivity	RNF111 - ring finger protein 111,
Immunogen	RNF111 (AAH60862.1, 1 a.a. - 985 a.a.) full-length human protein. MSQWTPEYKELYTLKVDKSEIPSDAPKTQESLKGILLHPEPIGAAKSFPAGVE MINSKVGNEFSHLCDSDSQKQEKEMNGNQEQEKSLVVRKKRKSQQAGPSYV QNCVKENQGILGLRQHLGTPSDEDNDSSFSDDLSSPSSSLHFGDSDTVTSD KEVSVRHSQTILNAKSRSHSARSHKWPRTEETESVSGLLMKRPCLHGSSLRRLP CRKRFVKNNSSQRTQKQKERILMQRKKREVLARRKYALLPSSSSSSSENDLSSE SSSSSSTEGEEDLFVSAENHQNNPAVPSGSDDEDVVIEASSTPQVTANEEIN VTSTDSEVEIVTVGESYRSRSTLGHRSRSHWSQGSSSSHASRPQEPNRNRSTV IQPLRQNAAEVVDLTVDEDEPTTVPTTSARMESQATSASINNSNPSTSEQASDT ASAVTSSQPSTVSETSATLTSTSTGTSIGDDSRRTTSSAVTETGPPAMPRLPS CCPQHSPCGGSSQNHHLGHPHTSCFQQHGHFQHGHGHGHHTHPHAPVPVSP SFSDPACPVERPPQVQAPCGANSSSGTSYHEQALPVDLSNSGIRSHGSGSFH GASAFDPCCPVSSSRAAIFGHQAAAAAPSQPLSSIDGYGSSMVAQPQPQPPP QPSLSSCRHYMPPPYASLTRPLHHQASACPHSHGNPPPQTQPPPQVDYVIPH PVHAFHSQISSHATSHPVAPPPPTHLASTAAPIQHLPPTHQPISHHIPATAPPA QRLHPHEVMQRMVQRRRMMQHPTRAHERPPPHPHRMHPNYGHGHHIHPV QTMSSHPRQAPERSAWELGIEAGVTAATYTPGALHPLAHYHAPPRLLHLQL GALPLMVPDMAGYPHIRYISSGLDGTSTFRGPFRGNFEELIHLEERLGNVNRGAS QGTIERCTYPHKYKKRKLHCKQDGEEGTEEDTEEKCTICLSILEEGEDVRRLLPC MHLFHQVCVDQWLITNKKCPICRVDIEAQLPSES
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500, Immunohistochemistry, Immunohistochemistry-Paraffin




Application Notes


Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control. Has also been used for immunohistochemistry (paraffin).

Images


Western Blot: RNF111 Antibody [H00054778-B01P] - Analysis of RNF111 expression in rat brain.



Immunohistochemistry-Paraffin: RNF111 Antibody [H00054778-B01P] - Analysis of purified antibody to RNF111 on formalin-fixed paraffin-embedded human prostate. (antibody concentration 3 ug/ml)



Western Blot: RNF111 Antibody [H00054778-B01P] - Analysis of RNF111 expression in transfected 293T cell line by RNF111 polyclonal antibody. Lane 1: RNF111 transfected lysate(108.35 KDa). Lane 2: Non-transfected lysate.

**Publications**

Zhou F, Drabsch Y, Dekker TJ et al. Nuclear receptor NR4A1 promotes breast cancer invasion and metastasis by activating TGF- β signalling. Nat Commun. 2014-03-03 [PMID: 24584437]

Sun Huaiyu, Liu Yijing, Hunter Tony. Multiple Arkadia/RNF111 structures coordinate its Polycomb body association and transcriptional control. Mol Cell Biol. 2014-08-01 [PMID: 24912682]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00054778-B01P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
NBP1-94140PEP	RNF111 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00054778-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

