

Product Datasheet

Asporin Antibody - Azide and BSA Free H00054829-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00054829-D01P

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00054829-D01P



H00054829-D01P

Asporin Antibody - Azide and BSA Free

Product Information

Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)

Product Description


Host	Rabbit
Gene ID	54829
Gene Symbol	ASPN
Species	Human, Mouse
Specificity/Sensitivity	Reacts with asporin.
Immunogen	ASPN (AAH63114.1, 1 a.a. - 384 a.a.) full-length human protein. MKEYVLLLFLALCSAKPFFSPSHIALKNMMLKDMEDTDDDDDDDDDDDDDDDDDD DEEDNSLFPTREPRSHFFPFDLFPMCPFGCQCYSRVVHCSDLGLTSVPTNIPF DTRMLDLQNNKIKEIKENDFKGLTSLYGLILNNNKLTKIHPKAFLTTKKLRRRLYS HNQLSEIPLNLPKSLAELRIHENKVKKIQKDTFKGMNALHVLEMSANPLDNNGIE PGA FEGVTVFHIRIAEAKLTSVPKGLPPTLLELHLDYNKISTVELEDFKRYKELQR LGLGNNKITDIENGLANIPRVREIHLENNKLKKIPSGLPELKYLQIIFLHSNSIARV GVNDFCPTVPMKKSLYSAISLFNNPVKYWEMQPATFRCVLSRMSVQLGNFG M
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details


Applications	Western Blot
Recommended Dilutions	Western Blot 1:100-1:2000
Application Notes	This antibody is reactive against tissue and transfected lysate in western blot, and as a detection antibody in ELISA.

Images


Western Blot: Asporin Antibody [H00054829-D01P] - Analysis of ASPN expression in mouse testis.




Western Blot: Asporin Antibody [H00054829-D01P] - Analysis of ASPN expression in HeLa.




Western Blot: Asporin Antibody [H00054829-D01P] - Analysis of ASPN expression in human kidney.



Western Blot: Asporin Antibody [H00054829-D01P] - Analysis of ASPN expression in human pancreas.



Western Blot: Asporin Antibody [H00054829-D01P] - Analysis of ASPN expression in transfected 293T cell line by ASPN polyclonal antibody. Lane 1: ASPN transfected lysate(43.90 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00054829-D01P

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NB100-1514PEP	Asporin Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00054829-D01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

