

# Product Datasheet

## Zinc finger protein 287 Antibody (3A2) - Azide and BSA Free H00057336-M08

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00057336-M08](http://www.novusbio.com/H00057336-M08)

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00057336-M08](http://www.novusbio.com/reviews/destination/H00057336-M08)




**H00057336-M08****Zinc finger protein 287 Antibody (3A2) - Azide and BSA Free**

| <b>Product Information</b>         |                                                                                                                                                                                                                                             |
|------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 mg                                                                                                                                                                                                                                      |
| <b>Concentration</b>               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.                                                                                                                            |
| <b>Storage</b>                     | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                                                                                                                                                |
| <b>Clonality</b>                   | Monoclonal                                                                                                                                                                                                                                  |
| <b>Clone</b>                       | 3A2                                                                                                                                                                                                                                         |
| <b>Preservative</b>                | No Preservative                                                                                                                                                                                                                             |
| <b>Isotype</b>                     | IgG2a Kappa                                                                                                                                                                                                                                 |
| <b>Purity</b>                      | IgG purified                                                                                                                                                                                                                                |
| <b>Buffer</b>                      | In 1x PBS, pH 7.4                                                                                                                                                                                                                           |
| <b>Product Description</b>         |                                                                                                                                                                                                                                             |
| <b>Description</b>                 | Novus Biologicals Mouse Zinc finger protein 287 Antibody (3A2) - Azide and BSA Free (H00057336-M08) is a monoclonal antibody validated for use in WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| <b>Host</b>                        | Mouse                                                                                                                                                                                                                                       |
| <b>Gene ID</b>                     | 57336                                                                                                                                                                                                                                       |
| <b>Gene Symbol</b>                 | ZNF287                                                                                                                                                                                                                                      |
| <b>Species</b>                     | Human                                                                                                                                                                                                                                       |
| <b>Specificity/Sensitivity</b>     | ZNF287 - zinc finger protein 287 (3A2)                                                                                                                                                                                                      |
| <b>Immunogen</b>                   | ZNF287 (NP_065704.1, 231 a.a. ~ 330 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.<br>PEWETKAQACTPVEDMSKLTKEETHTIKLEDSYDYDDRLERRGKGGFWKIHTD<br>ERGFSLKSVLSQEYDPTTEECLSKYDIYRNNFEKHSNLIVQFDTQLDN          |
| <b>Notes</b>                       | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                          |
| <b>Product Application Details</b> |                                                                                                                                                                                                                                             |
| <b>Applications</b>                | Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Sandwich ELISA                                                                                                                                                                |
| <b>Recommended Dilutions</b>       | Western Blot 1:500, ELISA, Immunocytochemistry/ Immunofluorescence, Sandwich ELISA                                                                                                                                                          |
| <b>Application Notes</b>           | Antibody Reactive Against Recombinant Protein with GST tag on ELISA and Western Blot. GST tag alone is used as a negative control.                                                                                                          |




## Images


Immunocytochemistry/Immunofluorescence: Zinc finger protein 287 Antibody (3A2) [H00057336-M08] - Immunofluorescence of monoclonal antibody to Zinc finger protein 287 on HeLa cell. [antibody concentration 40 ug/ml]



Immunocytochemistry/Immunofluorescence: Zinc finger protein 287 Antibody (3A2) [H00057336-M08] - Analysis of monoclonal antibody to ZNF287 on HeLa cell. Antibody concentration 40 ug/ml



Sandwich ELISA: Zinc finger protein 287 Antibody (3A2) [H00057336-M08] - Detection limit for recombinant GST tagged ZNF287 is 0.1 ng/ml as a capture antibody.



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to H00057336-M08**

---

|                    |                                                                |
|--------------------|----------------------------------------------------------------|
| HAF007             | Goat anti-Mouse IgG Secondary Antibody [HRP]                   |
| NB720-B            | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]        |
| NBP1-96981-0.5mg   | Mouse IgG2a Kappa Isotype Control (M2AK)                       |
| H00057336-P01-10ug | Recombinant Human Zinc finger protein 287 GST (N-Term) Protein |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00057336-M08](http://www.novusbio.com/reviews/submit/H00057336-M08)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

