

Product Datasheet

PCDHB16 Antibody (3H1) - Azide and BSA Free H00057717-M03

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00057717-M03

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00057717-M03



H00057717-M03

PCDHB16 Antibody (3H1) - Azide and BSA Free

Product Information	
Unit Size	0.05 mg
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3H1
Preservative	No Preservative
Isotype	IgG2b Kappa
Purity	Ascites
Buffer	In 1x PBS, pH 7.4


Product Description	
Description	Novus Biologicals Mouse PCDHB16 Antibody (3H1) - Azide and BSA Free (H00057717-M03) is a monoclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	57717
Gene Symbol	PCDHB16
Species	Human
Specificity/Sensitivity	PCDHB16 - protocadherin beta 16 (3H1)
Immunogen	PCDHB16 (AAH36062.1, 1 a.a. ~ 776 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MEIGWMHNRQRQVLVFFVLLSLSGAGAE LGSYSVVEETERGSFVANLGKDL GLGLTEMSTRKARIISQGNKQHLQLKAQTGDLLINEKLDREELCGPTEPCILHFQ VLMENPLEIFQAELRVIDINDHSPMFTEKEMILKIPENSPLGTEFPLNHALDLVDG SNNVQNYKISPSSHFRVLIHEFRDGRKYPELVLDKELDREEEPQLRLTLTALDG GSPPRSGTAQVRIEVDINDNAPEFEQPIYKVQIPENSPLGSLVATVSARDLDG GANGKISYTLFQPSSEDISKTLEVNPMTG EVRLRKQVDFEMVTSYEVRIKATDGG GLSGKCTLLLQVVDVNDNPPQVTMSALTSP IAENSPEIVVAVFSVSDPD SGNN GKTISSIQEDLPFLKPSVKNFYTLVTERALDREARAEYNITLTVTDMGTPRLKTE HNITVQISDVNDNAPTFTQTSYTLFVRENNSPALHIGSVSATDRD SGTNAQVTY SLLPPQDPHLPLASLV SINADNGHLFALRSLDY EALQAFEF RVGATDRGSPALS REALVRVLVLDANDNSPLVLYPLQNGSAPCTELVPRAAEPGYLVTKVVAVDGD SGQNAWLSYQLLKATEPGLFGVWAHNGEVRTARLLSERDAAKQRLVVLVKDN GEPPRSATATLHVLLVDGFSQPFLPLPEAAPGQTQANSLTVYLVVALASVSSLF LFSVLLFVAVRLCRRSRAASVGRCSMPEGPFPGR LVDVSGTGTLSQSYQYEV CLTGGSETSEFKFLKPIIPNFSP
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:500, ELISA
Application Notes	Antibody Reactive against Recombinant Protein with GST tag on ELISA and Western Blot. GST tag alone is used as a negative control.



Images

Western Blot: PCDHB16 Antibody (3H1) [H00057717-M03] - Detection against immunogen (111.1 kDa).



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00057717-M03

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)
DPI00	SLPI [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00057717-M03

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

