

Product Datasheet

PLEKHG2 Antibody - Azide and BSA Free H00064857-B01P-50ug

Unit Size: 50 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00064857-B01P-50ug

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00064857-B01P-50ug



H00064857-B01P-50ug

PLEKHG2 Antibody - Azide and BSA Free

Product Information	
Unit Size	50 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Novus Biologicals Mouse PLEKHG2 Antibody - Azide and BSA Free (H00064857-B01P-50ug) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	64857
Gene Symbol	PLEKHG2
Species	Human
Immunogen	PLEKHG2 (AAH13426.1, 1 a.a. - 896 a.a.) full-length human protein. MFPQNAKPGFKHAGSEGELYPPESQPPVSGSAPPEDLEDAGPPTLDPSGTSIT EEILELLNQRLRDPGPSTHDIPKFPGDSQVPGDSETLTFQALPSRDSSEEEEE EEEGLEMDERGPSPLHVLEGLESSIAAEMPSIPCLTKIPDVPNLPEIPSHCEIPE GSRLPSLSDISDVFEMPCLPAIPSVPNTPSLSSTPTLSCDSWLQGQLQEPAEAP ATRRELFSGSNPGKLGEPPSGGKAGPEEDEEGVSFTDFQPQDVTQHQGFPD ELAFRSCSEIRSAWQALEQQQLARPGFPEPLLILEDSDLGGDSGSGKAGAPSS ERTASRVRELARLYSERIQQMORAETRASANAPRRRPRVLAQPQPSPCLPQE QAEPGLLPAGFHVLCALAFPLTCAQESVPLGPAVWVQAAIPLSKQGGSPDGQ GLHVSNLPKQDLPGIHVSAATLLPEQGGSRHVQAPAATPLPKQEGPLHLQVPA LTTFSDQGHPEIQVPATTPLPEHKSHMVIPAPSTAFCEQGHCAIHVPTTPAL PKEICSDFTVSVTTPVVPKQEGHLDSESPTNIPLTKQGGSRDVQGPDPVCSQPIQ PLSWHGSSLDPQGGDTLPPLPCHLPDLQIPGTSPLPAHGSHLDHRIPANAPL SLSQELPDTQVPATTPLPLPQVLTDIWWQALPTSPKQGSPLDIQGPAAPPLPE PSLTDQVQKLTPSLEQKSLIDAHVPAATPLPERGGSLDIQGLSPTPVQTTMVL SKPGGSLASHVARLESSDLTPPHSPPPSSRQLLGPNAALSRYLAAASYISQSLA RRQGPGGGAPAAASRGSWSSAPTSSASSPPPPQPQPPPPAARRLSYATTVNIHV GGGRLRPAKAQVRLNHPALLASTQESMGLHRAQGAPDAPFHM
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot



Images

Western Blot: PLEKHG2 Antibody [H00064857-B01P-50ug] - analysis of PLEKHG2 expression in transfected 293T cell line by PLEKHG2 MaxPab polyclonal antibody. Lane 1: PLEKHG2 transfected lysate(98.56 KDa). Lane 2: Non-transfected lysate.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00064857-B01P-50ug

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
H00064857-P01-10ug	Recombinant Human PLEKHG2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00064857-B01P-50ug

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

