

Product Datasheet

ACD Antibody (1D8-1B6) - Azide and BSA Free H00065057-M02

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 8


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00065057-M02

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00065057-M02



H00065057-M02


ACD Antibody (1D8-1B6) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1D8-1B6
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Host	Mouse
Gene ID	65057
Gene Symbol	ACD
Species	Human
Reactivity Notes	Other species not tested.
Specificity/Sensitivity	ACD - adrenocortical dysplasia homolog (mouse)
Immunogen	ACD (AAH16904, 1 a.a. ~ 544 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MPGRCQSDAAMRVNGPASRAPAGWTSGLHTGPRAGRPRQAARGVVRGRGL LLRPRPAKELPLPRKGGAWAPAGNPGPLHPLGVAVGMAGSGRLVLRPWIRELI LGSETPSSPRAGQLLEVLQDAEAAVAGPSHAPDTSVVGATLLVSDGTHSVRCL VTREALDTSDWEEKEFGFRGTEGRLLLLQDCGVHVQVAEGGAPAEFYLVQVDR FSLLPTEQPRLRVPGCNQDLVDVQKKLYDCLEEHLSSESTSSNAGLSLSQLLDEM REDQEHQGALVCLAESCLTLEGPCAPPVTHWAASRCKATGEAVYTVPSSML CISENDQLILSSLGPCQRTQGPELPPDPALQDLSLTIASPPSSPSSSGTPALP GHMSSEESGTSISLLPALSLAAPDPGQRSSSQPSAICSAPATLTTPRSPHASRT PSSPLQSCTPSLSPRSHVPSPHQALVTRPQKPSLEFKEFVGLPCKNRPPFPRT GATRGAQEPCSVWEPPKRRHRDGSFAFYEYEPCTSLCARVQAARLPPQLMA WALHFLMDAQPGEPTPM
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP)
Recommended Dilutions	Western Blot 1:500, ELISA 1:100-1:2000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence 1:10-1:500, Immunoprecipitation 1:10-1:500, Immunohistochemistry-Paraffin 1:10-1:500, Chromatin Immunoprecipitation (ChIP) 1:10-1:500
Application Notes	Antibody reactive against cell lysate and recombinant protein for Western Blot. Has also been used for immunofluorescence, immunohistochemistry (paraffin), and ELISA. Use in chromatin immunoprecipitation reported in scientific literature (PMID 25233904)




Images


Western Blot: ACD Antibody (1D8-1B6) [H00065057-M02] - Western Blot analysis of ACD expression in IMR-32 (Cat # L008V1).




Immunocytochemistry/Immunofluorescence: ACD Antibody (1D8-1B6) [H00065057-M02] - Analysis of monoclonal antibody to ACD on HeLa cell. Antibody concentration 10 ug/ml.




Immunohistochemistry-Paraffin: ACD Antibody (1D8-1B6) [H00065057-M02] - Analysis of monoclonal antibody to ACD on formalin-fixed paraffin-embedded human heart tissue. Antibody concentration 2 ug/ml.



ELISA: ACD Antibody (1D8-1B6) [H00065057-M02] - Detection limit for recombinant GST tagged ACD is 0.3 ng/ml as a capture antibody.



Immunoprecipitation: ACD Antibody (1D8-1B6) [H00065057-M02] - Analysis of ACD transfected lysate using anti-ACD monoclonal antibody and Protein A Magnetic Bead, and immunoblotted with ACD MaxPab rabbit polyclonal antibody.



Publications

Sandhu R, Wei D, Sharma M et al. An N-terminal Flag-tag impairs TPP1 regulation of telomerase function *Biochem Biophys Res Commun* 2019-03-20 [PMID: 30885434]

Zanini IM, Sonesson C, Lorenzi LE, Azzalin CM. Human Cactin interacts with DHX8 and SRRM2 to assure efficient pre-mRNA splicing and sister chromatid cohesion. *J Cell Sci* 2017-02-15 [PMID: 28062851]

Frank AK, Tran DC, Qu RW et al. The Shelterin TIN2 Subunit Mediates Recruitment of Telomerase to Telomeres. *PLoS Genet.* 2015-07-31 [PMID: 26230315]

Nandakumar J, Bell CF, Weidenfeld I et al. The TEL patch of telomere protein TPP1 mediates telomerase recruitment and processivity. *Nature.* 2012-12-13 [PMID: 23103865] (Human)

Kocak H, Ballew BJ, Bisht K et al. Hoyerlaal-Hreidarsson syndrome caused by a germline mutation in the TEL patch of the telomere protein TPP1. *Genes Dev.* 2014-10-01 [PMID: 25233904] (Chemotaxis, Human)

Lakamp AS, Ouellette MM. A ssDNA Aptamer That Blocks the Function of the Anti-FLAG M2 Antibody. *J Nucleic Acids*;2011:720798. 2011-01-01 [PMID: 22013507]

Abreu E, Aritonovska E, Reichenbach P et al. TIN2-tethered TPP1 recruits human telomerase to telomeres in vivo. *Mol Cell Biol.* 2010-06-01 [PMID: 20404094]

Choi KH, Farrell AS, Lakamp AS et al. Characterization of the DNA binding specificity of Shelterin complexes. *Nucleic Acids Res.* 2011-11-01 [PMID: 21852327]

Details:
Tin2 (IMG-282) & POT1 (IMG-5343A). EMSA: HeLa cell line transfected with or without Flag-Tin2, Fig 4C.

**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00065057-M02

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
H00065057-P01-10ug	Recombinant Human ACD GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00065057-M02

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

