

Product Datasheet

TUBAL3 Antibody - Azide and BSA Free H00079861-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00079861-D01P

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00079861-D01P



H00079861-D01P

TUBAL3 Antibody - Azide and BSA Free

| Product Information | |
|---------------------|--|
| Unit Size | 0.1 mg |
| Concentration | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | No Preservative |
| Isotype | IgG |
| Purity | Protein A purified |
| Buffer | PBS (pH 7.4) |


| Product Description | |
|-------------------------|--|
| Description | Novus Biologicals Rabbit TUBAL3 Antibody - Azide and BSA Free (H00079861-D01P) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Rabbit |
| Gene ID | 79861 |
| Gene Symbol | TUBAL3 |
| Species | Human |
| Specificity/Sensitivity | Reacts with tubulin, alpha-like 3. |
| Immunogen | TUBAL3 (AAH98247.1, 1 a.a. - 406 a.a.) full-length human protein. MDQLENAMKMEHTNASFDTHFCETRAGKHVPRALFVDLEPTVIDGIRTGQHRSL FHPEQLLSGKEDAANNYARGRYSGSEVIDLVLERTRKLAEQCGGLQGFLIFR SFGGGTGSFGFTSLLMERLTGEYSRRTKLEFSVYPAPRISTAVVEPYNSVLTTSS TTEHTDCTFMVDNEAVYDICHRLKLGVECPHASINRLVVQVVSSITASLRFEGP LNVDLIEFQTNLVPYPRIHFPMTAFAPIVSADKAYHEQFSVSDITTACFESSNQL VKCDPRLGKYMACCLLYRGDVVPKEVNAAIAATKSRHSVQFVDWCPTGFKVGI NNRPPTVMPGGDLAKVHRSICMLSNTTAIVEAWARLDHKFDLMYAKRAFLHWY LREGMEEAEFLEAREDLAALERDYEEVAQSF |
| Notes | This product is produced by and distributed for Abnova, a company based in Taiwan. |

| Product Application Details | |
|-----------------------------|---|
| Applications | Western Blot |
| Recommended Dilutions | Western Blot |
| Application Notes | This antibody is reactive against transfected lysate in western blot, and as a detection antibody in ELISA. |



Images

Western Blot: TUBAL3 Antibody [H00079861-D01P] - Analysis of TUBAL3 expression in transfected 293T cell line by TUBAL3 polyclonal antibody. Lane 1: TUBAL3 transfected lysate(45.50 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00079861-D01P

| | |
|--------------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| H00079861-P01-10ug | Recombinant Human TUBAL3 GST (N-Term) Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00079861-D01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

