

Product Datasheet

IFT20 Antibody (3F3) - Azide and BSA Free H00090410-M02

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00090410-M02

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00090410-M02



H00090410-M02


IFT20 Antibody (3F3) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3F3
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Novus Biologicals Mouse IFT20 Antibody (3F3) - Azide and BSA Free (H00090410-M02) is a monoclonal antibody validated for use in WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	90410
Gene Symbol	IFT20
Species	Human
Immunogen	IFT20 (NP_777547.1, 1 a.a. ~ 148 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MAKDILGEAGLHFDELNKLRLDPEVTQQTIELKEECKDFVDKIGQFQKIVGGI ELVDQLAKEAENEKMKSLAVSPRLECTGAISAHCKLCLSDSSDSPTSPSRVGG TTGHRCSLAQIYSKAERSSTAATSSPNSRKENAARKVSG
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Sandwich ELISA
Recommended Dilutions	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Sandwich ELISA
Application Notes	This antibody is useful for ELISA, Western Blot




Images


Western Blot: IFT20 Antibody (3F3) [H00090410-M02] - Analysis of IFT20 expression in transfected 293T cell line by IFT20 monoclonal antibody (M02), clone 3F3. Lane 1: IFT20 transfected lysate(16 KDa). Lane 2: Non-transfected lysate.



Immunocytochemistry/Immunofluorescence: IFT20 Antibody (3F3) [H00090410-M02] - Analysis of monoclonal antibody to IFT20 on HeLa cell . Antibody concentration 10 ug/ml.



Sandwich ELISA: IFT20 Antibody (3F3) [H00090410-M02] - Detection limit for recombinant GST tagged IFT20 is 0.1 ng/ml as a capture antibody.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00090410-M02

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
NBP2-23152	Recombinant Human IFT20 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00090410-M02

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

