

# Product Datasheet

## Argonaute 3 Antibody - Azide and BSA Free H00192669-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00192669-D01P](http://www.novusbio.com/H00192669-D01P)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00192669-D01P](http://www.novusbio.com/reviews/destination/H00192669-D01P)



**H00192669-D01P**

Argonaute 3 Antibody - Azide and BSA Free

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 mg                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                     |
| Clonality           | Polyclonal                                                                                                       |
| Preservative        | No Preservative                                                                                                  |
| Isotype             | IgG                                                                                                              |
| Purity              | Protein A purified                                                                                               |
| Buffer              | PBS (pH 7.4)                                                                                                     |


| Product Description     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description             | Novus Biologicals Rabbit Argonaute 3 Antibody - Azide and BSA Free (H00192669-D01P) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Host                    | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Gene ID                 | 192669                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Gene Symbol             | AGO3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Species                 | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Specificity/Sensitivity | Rabbit Polyclonal to EIF2C3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Immunogen               | EIF2C3 (NP_803171.1, 1 a.a. - 626 a.a.) full-length human protein.<br>MCEVLDIHÑIDEQPRPLTDSHRVKFTKEIKGLKVEVTHCGTMRRKYRVCNVTRR<br>PASHQTFPLQLENGQTVERTVAQYFREKYTLQLKYPHLPCLQVGQEQKHTYLP<br>LEVCNIVAGQRCIKKLTDNQTSTMIKATARSAPDRQEEISRLVRSANYETDPFV<br>QEFQFKVRDEMAHVTGRVLPAPMLQYGGRNRTVATPSHGVWDMRGKQFHT<br>GVEIKMWAIACFATQRQCREEILKGFTDQLRKISKDAGMPIQGQPCFCKYAQG<br>ADSVEPMFRHLKNTYSGLQLIIVILPGKTPVYAEVKRVGDTLLGMATQCVQVK<br>VIKTSPQTLNLCLKINVKLGGINNILVPHQRPSVFQQPVIFLGADVTHPPAGDG<br>KKPSIAAVVGSMDAHPSTRYCATVRVQRPRQEIIQDLASMVRELLIQFYKSTRFK<br>PTRIIFYRDGVSEGQFRQVLYYELLAIREACISLEKDYQPGITYIVVQKRHHTRLF<br>CADRTERVGRSGNIPAGTTVDTDITHPYEFDYLCSHAGIQGTSRPSHYHVLW<br>DDNCFTADELQLLTYQLCHTYVRCTRSVSIPAPAYY AHLVAFRARYHLVDKEHD<br>SAEGSHVSGQSNGRDPQALAKAVQIHQDTRLRTMYFA |
| Notes                   | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

| Product Application Details |                                                                                                                                                                                   |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Western Blot                                                                                                                                                                      |
| Recommended Dilutions       | Western Blot 1:500                                                                                                                                                                |
| Application Notes           | Antibody reactive against recombinant protein with GST tag on ELISA and western blot and also on transfected lysate in western blot. GST tag alone is used as a negative control. |



## Images

Western Blot: Argonaute 3 Antibody [H00192669-D01P] - Analysis of EIF2C3 expression in transfected 293T cell line by EIF2C3 polyclonal antibody. Lane 1: EIF2C3 transfected lysate (71.20 KDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00192669-D01P**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP2-57983PEP | Argonaute 3 Recombinant Protein Antigen             |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00192669-D01P](http://www.novusbio.com/reviews/submit/H00192669-D01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

