

Product Datasheet

PI4KAP2 Antibody - Azide and BSA Free H00375133-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00375133-D01P

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00375133-D01P



H00375133-D01P

PI4KAP2 Antibody - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Novus Biologicals Rabbit PI4KAP2 Antibody - Azide and BSA Free (H00375133-D01P) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	375133
Gene Symbol	PI4KAP2
Species	Human
Specificity/Sensitivity	Reacts with phosphatidylinositol 4-kinase, catalytic, alpha pseudogene 2.
Immunogen	PI4KAP2 (AAI41653.1, 1 a.a. - 592 a.a.) full-length human protein. MTVEQKFGFLFAEIKEADPLAASEASQPKPCPEVTPHYIWIDFLVQRFEIAKYC SSDQVEIFSSLLQRSMSLNIGRAKGSMMNRHVAAIGPRFKLLTLGLSLLHADVVP NATIRNVLREKIYSTAFDYFSCPPKFPTQGEKRLREDISIMIKFWTAMFSDKKYL TASQLVPPADIGDLEQLVEENTGSLSGPAKDFYQREFDFFNKITNVSIIKPYP KGDERKKACLSALSEVTVQPGCSLPSNPEAIVLDVDYKSGTPMQSAAKAPYLA KFKVKRCGVSELEKEGLRCRSDSEDECSTQEADGQKISWQAAIFKLGDDCRQ KSYWGARTPTDRILRLPASQDMLALQIIDLFKNIFQLVGLDLFVFPYRVVATAPG CGVIECIPDCTSRDQLGRQTDGFMYYDYFTRQYGDDESTLAFQQACYNFIRSM YSLLLFLLQIKDRHNGNIMLDKKGHIIHIDFGFMFESSPGGNLGWEPDIKLTDEM VMIMGGKMEATPFKWFMEMCVQATWLCGEPGPYMDVVVSLVTIMLDTGLPCF RGQTIKFLKHRFSPNMTEREANFIMKVIQSCFLSNRSRTYNMIQYYQNDIPY
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	This antibody is reactive against tissue and transfected lysate in western blot, and as a detection antibody in ELISA.




Images

Western Blot: PI4KAP2 Antibody [H00375133-D01P] - Analysis of PI4KAP2 expression in human liver.



Western Blot: PI4KAP2 Antibody [H00375133-D01P] - Analysis of PI4KAP2 expression in transfected 293T cell line by PI4KAP2 polyclonal antibody. Lane 1: PI4KAP2 transfected lysate (65.12 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00375133-D01P

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00375133-P02-10ug	Recombinant Human PI4KAP2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00375133-D01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

