

Product Datasheet

CATIP Antibody (2C3) - Azide and BSA Free H00375307-M01

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00375307-M01

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00375307-M01




H00375307-M01

CATIP Antibody (2C3) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2C3
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Host	Mouse
Gene ID	375307
Gene Symbol	CATIP
Species	Human
Reactivity Notes	Human. Other species not tested.
Specificity/Sensitivity	MGC50811 - hypothetical gene supported by BC052750
Immunogen	MGC50811 (AAH52750, 1 a.a. ~ 387 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MSSKVYSTGSRKDHQPSGPECLPLPEANAEIDFLSSLHKEELQMLFFSETLA MVSDTGEPQGELTIEVQRGKYQEKLGMALTYCLFVHASSRGFLDKMLCGNSLL GYLSEKLELMEQHSQDFIKFLILPMEKMSLLKQDDQLAVTRSIKEGEEVKTGV TSFPWSSIKGFISEAANLVLLRVMWRRMVPSNARFLTLDTEGKLCYLTQNLG FQTIQVDHQAEVFIVEQTVHAEEGIPMSCQYLLSDGHLAKRIQVVGSPGCCIIIT KMPILREEDIEPRPVFEKKPLVWEEDMELYSKFLDRKEELRLGHASYLRQHPE AHALISDFLLFLLLRQPEDVVTFAAEFFGPFDPWRPSSPALGSSHRPNPFRSLE PEGDARSGAA
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	ELISA, Sandwich ELISA
Recommended Dilutions	ELISA, Sandwich ELISA
Application Notes	This antibody is useful for ELISA.

Images

Sandwich ELISA: C2orf62 Antibody (2C3) [H00375307-M01] - Detection limit for recombinant GST tagged C2orf62 is 0.3 ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00375307-M01

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
H00375307-P01-10ug	Recombinant Human CATIP GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00375307-M01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

