

Product Datasheet

DNMT3A Assay Kit (Colorimetric) KA0683

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/KA0683

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/KA0683



KA0683**DNMT3A Assay Kit (Colorimetric)**

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.
Product Description	
Description	DNMT3A Assay Kit is use for measuring the amount of DNMT3A.
Gene ID	1788
Gene Symbol	DNMT3A
Kit Components	DA1 (10X Wash Buffer), DA2 (Assay Buffer), DA3 (DNMT3A Standard, 20 µg/mL), DA4 (Capture Antibody, 100 µg/mL), DA5 (Detection Antibody 200 µg/mL), DA6 (Developing Solution), DA7 (Stop Solution), BB Blocking Buffer, 8-Well Assay Strip (with Frame)
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Standard Curve Range	1 to 20 ug sample
Sensitivity	1 ug sample
Assay Type	Colorimetric
Suitable Sample Type	Nuclear extracts
Sample Volume	2- 6 uL
Product Application Details	
Applications	Quantification
Recommended Dilutions	Quantification



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to KA0683

H00001788-Q01-10ug	Recombinant Human DNMT3A GST (N-Term) Protein
NB100-56519	DNMT1 Antibody (60B1220.1) - BSA Free
NB120-13888	DNMT3A Antibody (64B1446) - BSA Free
AF1759	Oct-3/4 Antibody [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/KA0683

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

