

Product Datasheet

KAT1 Inhibitor Screening Assay Kit KA1318

Unit Size: 1 Kit

Store at -20C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/KA1318

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/KA1318




KA1318**KAT1 Inhibitor Screening Assay Kit**

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at -20C.
Product Description	
Description	HAT Inhibitor Screening Assay Kit is a fluorescence-based method for screening pCAF HAT inhibitors.
Gene ID	8520
Gene Symbol	HAT1
Kit Components	HAT Assay Buffer (5X), HAT Acetyl CoA, Histone Acetyltransferase (pCAF), HAT Peptide, HAT Stop Reagent, HAT Developer, 96-Well Plate (white), 96-Well Cover Sheet
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Suitable Sample Type	Compounds expected to inhibit PCAF
Sample Volume	5 uL/well. All samples should be assayed in duplicate at minimum.
Product Application Details	
Applications	Functional
Recommended Dilutions	Functional

Images

KAT1 Assay/Inhibitor Screening Assay Kit [KA1318] - Interference of garcinol in the assay.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to KA1318

NBP1-87835PEP	KAT1 Recombinant Protein Antigen
AF887	AKT [p Ser473] Antibody [Unconjugated] - Pan Specific
NBP1-87835	KAT1 Antibody
NB500-171	Histone H3 Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/KA1318

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

