

Product Datasheet

RSK1 [p Ser380] Antibody (1024F) - Azide and BSA Free MAB79672

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/MAB79672

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/MAB79672



MAB79672

RSK1 [p Ser380] Antibody (1024F) - Azide and BSA Free

Product Information	
Unit Size	100 ug
Concentration	1.0 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1024F
Preservative	No Preservative
Isotype	IgG
Purity	Protein A or G purified
Buffer	PBS


Product Description	
Host	Rabbit
Gene ID	6195
Gene Symbol	RPS6KA1
Species	Human
Immunogen	Phosphopeptide containing the human RSK1 S380 site Accession # Q15418

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 0.25 ug/ml, Immunocytochemistry/ Immunofluorescence




Images

Western Blot: RSK1 [p Ser380] Antibody (1024F) [MAB79672] - Analysis of HeLa +/- PMA cell lysates with 0.25 ug/ml of hRSK (S380) Monoclonal Antibody followed by HAF008. A specific band was detected for hRSK (S380) at approximately 90 kDa (as indicated). This experiment was conducted under reducing conditions and using Immunoblot Buffer Group 1.



Immunocytochemistry/Immunofluorescence: RSK1 [p Ser380] Antibody (1024F) [MAB79672] - RSK1 [pSer380] was detected in immersion fixed HeLa cells treated with 200 nm of PMA using RSK1 [pSer380] Monoclonal Antibody at 1 ug/ml dilution. Cells were stained using the NorthernLights 557-conjugated Anti-Rabbit IgG Secondary Antibody (red) and counterstained with DAPI (blue). Specific staining was localized to the cytoplasm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to MAB79672

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-89647PEP	RSK1 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/MAB79672

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

