

Product Datasheet

TMIGD2 Antibody (953730) - Azide and BSA Free MAB83161

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/MAB83161

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/MAB83161



MAB83161

TMIGD2 Antibody (953730) - Azide and BSA Free

Product Information	
Unit Size	100 ug
Concentration	1.0 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	953730
Preservative	No Preservative
Isotype	IgG2b
Purity	Protein A or G purified
Buffer	PBS


Product Description	
Host	Mouse
Gene ID	126259
Gene Symbol	TMIGD2
Species	Human
Immunogen	Partial recombinant human CD28H protein (amino acids 23-150) [UniProt Q96BF3]

Product Application Details	
Applications	Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Paraffin, CyTOF-ready, Western Blot (Negative)
Recommended Dilutions	Flow Cytometry, Immunohistochemistry 5 ug/ml, Immunohistochemistry-Paraffin, Western Blot (Negative), CyTOF-ready
Application Notes	This antibody is Cytof ready.




Images

Immunohistochemistry-Paraffin: TMIGD2 Antibody (953730) - Azide Free [MAB83161] - TMIGD2 was detected in immersion fixed paraffin-embedded sections of human thymus using tmigd2 Monoclonal Antibody at 5 ug/ml. Tissue was stained using Goat anti-Mouse IgG Secondary Antibody [HRP Polymer] (brown) and counterstained with hematoxylin (blue). Specific staining was localized to plasma membrane in cancer cells.



Flow Cytometry: TMIGD2 Antibody (953730) - Azide Free [MAB83161] - Detection of TMIGD2 in HEK293 Human Cell Line Transfected with Human TMIGD2 and eGFP by Flow Cytometry. HEK293 human embryonic kidney cell line transfected with either (A) human TMIGD2 or (B) irrelevant transfectants and eGFP was stained with Mouse Anti-Human TMIGD2 Monoclonal Antibody (Catalog # MAB83161) followed by Allophycocyanin-conjugated Anti-Mouse IgG Secondary Antibody (Catalog # F0101B). Quadrant markers were set based on control antibody staining (Catalog # MAB0041).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to MAB83161

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)
NBP1-81164PEP	TMIGD2 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/MAB83161

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

