

# Product Datasheet

## MLKL [p Thr357] Antibody (954724) [Allophycocyanin] MAB91871APC

Unit Size: 100 ul

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/MAB91871APC](http://www.novusbio.com/MAB91871APC)

Updated 10/30/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/MAB91871APC](http://www.novusbio.com/reviews/destination/MAB91871APC)



**MAB91871APC**

MLKL [p Thr357] Antibody (954724) [Allophycocyanin]

| Product Information         |                                                                                                                   |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                            |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                       |
| Storage                     | Store at 4C in the dark.                                                                                          |
| Clonality                   | Monoclonal                                                                                                        |
| Clone                       | 954724                                                                                                            |
| Preservative                | 0.05% Sodium Azide                                                                                                |
| Isotype                     | IgG2a                                                                                                             |
| Conjugate                   | Allophycocyanin                                                                                                   |
| Purity                      | Protein A or G purified                                                                                           |
| Buffer                      | PBS                                                                                                               |
| Product Description         |                                                                                                                   |
| Description                 | >90% by SDS PAGE                                                                                                  |
| Host                        | Mouse                                                                                                             |
| Gene ID                     | 197259                                                                                                            |
| Gene Symbol                 | MLKL                                                                                                              |
| Species                     | Human                                                                                                             |
| Immunogen                   | Synthetic peptide from human MLKL phosphorylated at T357 [UniProt Q8NB16]                                         |
| Product Application Details |                                                                                                                   |
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Flow Cytometry (Negative), Immunocytochemistry |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, Flow Cytometry (Negative), Immunocytochemistry |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                            |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to MAB91871APC**

|                    |                                                       |
|--------------------|-------------------------------------------------------|
| IC003A             | Mouse IgG2a Isotype Control (20102) [Allophycocyanin] |
| H00197259-Q01-10ug | Recombinant Human MLKL GST (N-Term) Protein           |
| 210-TA-005         | TNF-alpha [Unconjugated]                              |
| H00197259-T01      | MLKL 293T Cell Transient Overexpression Lysate        |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/MAB91871APC](http://www.novusbio.com/reviews/submit/MAB91871APC)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

