

Product Datasheet

DAP10/HCST Antibody (982101) [Unconjugated] - BSA Free MAB97861

Unit Size: 100 ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/MAB97861

Updated 2/21/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/MAB97861



MAB97861

DAP10/HCST Antibody (982101) [Unconjugated] - BSA Free


Product Information	
Unit Size	100 ug
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	982101
Preservative	0.02% Sodium Azide
Isotype	IgG2a
Conjugate	Unconjugated
Purity	Protein A or G purified
Buffer	PBS

Product Description	
Host	Mouse
Gene ID	10870
Gene Symbol	HCST
Species	Human
Immunogen	Synthetic peptide containing human DAP10/HCST

Product Application Details	
Applications	Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Flow Cytometry 1:10 - 1:1000, Immunohistochemistry 5 - 10 ug/ml, Immunohistochemistry-Paraffin 5 - 10 ug/ml

Images

Immunohistochemistry-Paraffin: DAP10/HCST Antibody (982101) [Unconjugated] [MAB97861] - DAP10/HCST was detected in immersion fixed paraffin-embedded sections of human kidney using Mouse Anti-Human DAP10/HCST Monoclonal Antibody (Catalog # MAB97861) at 5 ug/mL for 1 hour at room temperature followed by incubation with the Anti-Mouse IgG VisUCyte HRP Polymer Antibody (Catalog # VC001). Before incubation with the primary antibody, tissue was subjected to heat-induced epitope retrieval using Antigen Retrieval Reagent-Basic (Catalog # CTS013). Tissue was stained using DAB (brown) and counterstained with hematoxylin (blue). Specific staining was localized to cell nuclei.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to MAB97861

MAB003	Mouse IgG2a Isotype Control (20102) [Unconjugated]
NBP2-61459-50ug	Recombinant Human DAP10/HCST hIgG-His Protein
210-TA-005	TNF-alpha [Unconjugated]
NBL1-11475	DAP10/HCST Overexpression Lysate

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/MAB97861

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

