

Product Datasheet

DAP3 Antibody NB100-1032

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-1032

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-1032



NB100-1032**DAP3 Antibody**

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris saline (20 mM Tris pH 7.3, 150 mM NaCl), 0.5% BSA


Product Description	
Description	Novus Biologicals Goat DAP3 Antibody (NB100-1032) is a polyclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. Anti-DAP3 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Goat
Gene ID	7818
Gene Symbol	DAP3
Species	Human
Specificity/Sensitivity	This antibody is expected to recognize isoform 1 (NP_387506.1), isoform 2 (NP_001186779.1) and isoform 3 (NP_001186780.1). Reported variants represent identical protein (NP_387506.1; NP_004623.1; NP_001186778.1).
Immunogen	Peptide with sequence NPSLLERHCAYL corresponding to C-Terminus according to NP_387506.1, NP_001186779.1, NP_001186780.1.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Peptide ELISA
Recommended Dilutions	Western Blot 0.1 - 1 ug/mL, Immunohistochemistry, Immunocytochemistry/Immunofluorescence 10 ug/mL, Immunohistochemistry-Paraffin 5 ug/mL, Peptide ELISA Detection limit 1:64000
Application Notes	WB: Approx 40 kDa band observed in lysates of cell lines HeLa and HepG2 and in Human Kidney lysates (calculated MW of 41 kDa according to NP_001186779.1 and NP_001186780). Primary incubated 1 hour at room temperature.




Images


Western Blot: DAP3 Antibody [NB100-1032] - Staining (0.3 ug/mL) of Human Kidney lysate (RIPA buffer, 30 ug total protein per lane). Detected by chemiluminescence.




Immunocytochemistry/Immunofluorescence: DAP3 Antibody [NB100-1032] - Analysis of paraformaldehyde fixed MCF7 cells, permeabilized with 0.15% Triton. Primary incubation 1hr (10 ug/mL) followed by Alexa Fluor 488 secondary antibody (2 ug/mL), showing nuclear staining. The nuclear stain is DAPI (blue). Negative control: Unimmunized goat IgG (10 ug/mL) followed by Alexa Fluor 488 secondary antibody (2 ug/mL).




Western Blot: DAP3 Antibody [NB100-1032] - Staining of HeLa (A) and HepG2 (B) cell lysate (RIPA buffer, 30 ug total protein per lane). Antibody at 0.3 ug/mL. Detected by chemiluminescence.



Immunohistochemistry-Paraffin: DAP3 Antibody [NB100-1032] - Negative Control showing staining of paraffin embedded Human Kidney, with no primary antibody.



Immunohistochemistry-Paraffin: DAP3 Antibody [NB100-1032] - (5ug/ml) staining of paraffin embedded Human Kidney. Heat induced antigen retrieval with citrate buffer pH 6, HRP-staining.



Publications

Kissil JL, Deiss LP, Bayewitch M et al. Isolation of DAP3, a novel mediator of interferon-gamma-induced cell death. J Biol Chem 1995-11-17 [PMID: 7499268]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-1032

NB800-PC1	HeLa Whole Cell Lysate
HAF017	Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
HAF109	Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB410-28088-1mg	Goat IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-1032

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



