

# Product Datasheet

## ING5 Antibody [DyLight 650]

### NB100-1681C

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-1681C](http://www.novusbio.com/NB100-1681C)

Updated 9/21/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-1681C](http://www.novusbio.com/reviews/destination/NB100-1681C)



**NB100-1681C**

ING5 Antibody [DyLight 650]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                                                    |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Conjugate                   | DyLight 650                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Purity                      | Immunogen affinity purified                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Host                        | Goat                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Gene ID                     | 84289                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Gene Symbol                 | ING5                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Species                     | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Reactivity Notes            | Reactivity occurs against human p28 ING5 protein. The sequence is present in both isoforms 1 and 2. However, 100% homology is on record for this protein from human and chimpanzee. Cross reactivity with p28 ING5 protein from frog, mouse, rat and chicken may not occur as sequence homology varies by one or more amino acids residues in this sequence as indicated by BLAST analysis. Reactivity with p28 ING5 proteins from other sources is not known. |
| Specificity/Sensitivity     | Reactivity occurs against human p28 ING5 protein. The sequence is present in both isoforms 1 and 2.                                                                                                                                                                                                                                                                                                                                                            |
| Immunogen                   | This affinity purified ING5 Antibody was prepared from whole goat serum produced by repeated immunizations with a synthetic peptide corresponding to an internal region near aa 125-150 of Human p28 ING5 protein (Inhibitor of growth family, member 5). (Uniprot: Q8WYH8)                                                                                                                                                                                    |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                                                                                                                                              |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Applications                | Western Blot, ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Recommended Dilutions       | Western Blot, ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                                                         |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NB100-1681C**

---

|                    |                                             |
|--------------------|---------------------------------------------|
| NB410-28088C-500ul | Goat IgG Isotype Control [DyLight 650]      |
| H00084289-P01-10ug | Recombinant Human ING5 GST (N-Term) Protein |
| 2695-SE-010        | Airway Trypsin-like Protease/HAT/TMPRSS11D  |
| H00084289-Q01-10ug | Recombinant Human ING5 GST (N-Term) Protein |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-1681C](http://www.novusbio.com/reviews/submit/NB100-1681C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

