

Product Datasheet

DDX21 Antibody NB100-1717

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-1717

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-1717



NB100-1717

DDX21 Antibody


| Product Information | |
|-------------------------|-----------------------------|
| Unit Size | 100 ul |
| Concentration | 0.2 mg/ml |
| Storage | Store at 4C. Do not freeze. |
| Clonality | Polyclonal |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG |
| Purity | Immunogen affinity purified |
| Buffer | TBS and 0.1% BSA |
| Target Molecular Weight | 87 kDa |

| Product Description | |
|---------------------|---|
| Host | Rabbit |
| Gene ID | 9188 |
| Gene Symbol | DDX21 |
| Species | Human |
| Immunogen | The immunogen recognized by this antibody maps to a region between residues 50 and 100 of human DEAD/H (Asp-Glu-Ala-Asp/His) Box Polypeptide 21 using the numbering given in entry NP_004719.2 (GeneID 9188). |

| Product Application Details | |
|-----------------------------|---|
| Applications | Western Blot, Immunoprecipitation, Knockdown Validated |
| Recommended Dilutions | Western Blot 1:2000-1:10000, Immunoprecipitation 1-4 ug/mg of lysate, Knockdown Validated |

Images

Western Blot: DDX21 Antibody [NB100-1717] - Detection of Human DDX21 on HeLa whole cell lysate using NB100-1717. DDX21 was also immunoprecipitated using rabbit anti-DDX21 antibodies NB100-1716 and NB100-1718.

**Publications**

McRae EKS, Booy EP, Moya-Torres A et al. Human DDX21 binds and unwinds RNA guanine quadruplexes Nucleic Acids Res. 2017-06-19 [PMID: 28472472] (IP, KD, WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-1717

| | |
|------------|---|
| NBL1-09787 | DDX21 Overexpression Lysate |
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-1717

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

