

Product Datasheet

SMC3 Antibody - BSA Free

NB100-207

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 9

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-207

Updated 2/21/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-207



NB100-207

SMC3 Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)

Product Description

Host	Rabbit
Gene ID	9126
Gene Symbol	SMC3
Species	Human, Mouse, C. elegans
Reactivity Notes	Orangutan (100%). C. elegans reactivity reported in scientific literature (PMID: 30563860).
Immunogen	The immunogen recognized by this antibody maps to a region between residue 1175 and the C-terminus (residue 1217) of human Structural Maintenance of Chromosomes 3 using the numbering given in entry NP_005436.1 (GeneID 9126).


Product Application Details

Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Western Blot 1:10000 - 1:25000, Immunohistochemistry, Immunoprecipitation 2-10 ul/mg, Immunohistochemistry-Paraffin 1:100-1:500
Application Notes	Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.




Images


Western Blot: SMC3 Antibody [NB100-207] - Whole cell lysate (10 ug) from HEK293T, HeLa, Jurkat, LNCaP, U2OS, A-549, TCMK-1, and NIH 3T3 cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-SMC3 antibody used for WB at 0.04 ug/ml. Detection: Chemiluminescence with an exposure time of 3 seconds.




Immunohistochemistry-Paraffin: SMC3 Antibody [NB100-207] - Section of mouse renal cell carcinoma. Antibody: Affinity purified rabbit anti-SMC3 used at 1:5,000 (0.2ug/ml). Detection: DAB



Immunohistochemistry-Paraffin: SMC3 Antibody [NB100-207] - Section of human breast carcinoma. Antibody: Affinity purified rabbit anti-SMC3 used at 1:5,000 (0.2ug/ml). Detection: DAB



SMC3 Antibody [NB100-207] - Detection of human SMC3 by western blot of immunoprecipitates. Samples: Whole cell lysate (1.0 mg per IP reaction; 20% of IP loaded) from HeLa, cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-SMC3 antibody NB100-207 (lot NB100-207-9) used for IP at 3 ug per reaction. SMC3 was also immunoprecipitated by a previous lot of this antibody (lot NB100-207-8). For blotting immunoprecipitated SMC3, NB100-207 was used at 1 ug/ml. Detection: Chemiluminescence with an exposure time of 10 seconds.



Publications

Tao, H, Lambert, J P Et al. IRX3/5 regulate mitotic chromatid segregation and limb bud shape. Development 2020-10-05 [PMID: 32907847] (Chemotaxis, Human)

Bowman R, Balukof N, Ford T, Smolikove S. A Novel Role for a-importins and Akirin in Establishment of Meiotic Sister Chromatid Cohesion in *Caenorhabditis elegans* Genetics 2018-12-18 [PMID: 30563860] (IF/IHC, C. elegans)

Deardorff MA, Bando M, Nakato R et al. HDAC8 mutations in Cornelia de Lange syndrome affect the cohesin acetylation cycle. Nature 2012-09-01 [PMID: 22885700]

Shiomi Y, Hayashi A, Ishii T et al. Two different replication factor C proteins, Ctf18 and RFC1, separately control PCNA-CRL4Cdt2-mediated Cdt1 proteolysis during S phase and following UV irradiation. Mol Cell Biol. 2012-06-01 [PMID: 22493068]

Leman AR, Noguchi C, Lee CY et al. Human Timeless and Tipin stabilize replication forks and facilitate sister-chromatid cohesion. J Cell Sci 2010-03-01 [PMID: 20124417]

Terret ME, Sherwood R, Rahman S et al. Cohesin acetylation speeds the replication fork. Nature 2009-11-01 [PMID: 19907496]

Adelfalk C, Janschek J, Revenkova E et al. Cohesin SMC1beta protects telomeres in meiocytes. J Cell Biol 2009-10-01 [PMID: 19841137]

Rankin S, Ayad NG, Kirschner MW et al. Sororin, a substrate of the anaphase-promoting complex, is required for sister chromatid cohesion in vertebrates. Mol Cell 2005-04-01 [PMID: 15837422]

Revenkova E, Eijpe M, Heyting C et al. Cohesin SMC1 beta is required for meiotic chromosome dynamics, sister chromatid cohesion and DNA recombination. Nat Cell Biol 2004-06-01 [PMID: 15146193]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-207

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00009126-Q01-10ug	Recombinant Human SMC3 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-207

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

