

# Product Datasheet

## LC3B Antibody [mFluor Violet 450 SE] NB100-2220MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-2220MFV450](http://www.novusbio.com/NB100-2220MFV450)

Updated 9/20/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-2220MFV450](http://www.novusbio.com/reviews/destination/NB100-2220MFV450)



**NB100-2220MFV450**

LC3B Antibody [mFluor Violet 450 SE]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	mFluor Violet 450 SE
Purity	Immunogen affinity purified
Buffer	50mM Sodium Borate
Product Description	
Host	Rabbit
Gene ID	81631
Gene Symbol	MAP1LC3B
Species	Human, Mouse, Rat, Porcine, Alligator, Avian, Bacteria, Bovine, Canine, Chicken, Chinese Hamster, Guinea Pig, Hamster, Invertebrate, Monkey, Primate, Rabbit, Golden Syrian Hamster, Zebrafish
Reactivity Notes	Use in Rat reported in scientific literature (PMID:34622072). Use in Chinese Hamster reported in scientific literature (PMID:34332287). Mouse reactivity reported in scientific literature (PMID:32814898). Bovine reactivity reported in scientific literature (PMID: 24895572). Primate reactivity reported in scientific literature (PMID: 25142602). Canine reactivity reported in scientific literature (PMID: 25839646). Avian reactivity reported in scientific literature (PMID: 29546310). Hamster reactivity reported in scientific literature (PMID: 26423766). Rabbit reactivity reported in scientific literature (PMID: 26497211). The mouse detection has been reported to be weaker than the human. Immunogen displays the following percentage of sequence identity for non-tested species: Xenopus 84%. Invertebrate reactivity reported in scientific literature (PMID: 26716072 ). Monkey reactivity reported in scientific literature (PMID: 30324853). Guinea pig reactivity reported from a verified customer review. Rat reactivity reported in scientific literature (PMID: 30744518). Bacteria reactivity reported in scientific literature (PMID: 28783414). Chicken reactivity reported in scientific literature (PMID: 30649814). Porcine reactivity reported in scientific literature (PMID: 30789643). Golden Syrian hamster reactivity reported in scientific literature (PMID: 23180219). Zebrafish reactivity reported in scientific literature (PMID: 29185873). Human reactivity reported in scientific literature (PMID:33086573). Use in Alligator reported in scientific literature (PMID:32061056).
Marker	Autophagosome Marker
Immunogen	Polyclonal LC3B Antibody was made to a synthetic peptide made to an N-terminal portion of the human LC3B protein sequence (between residues 1-100). [UniProt# Q9GZQ8]
Notes	mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Product Application Details	



<b>Applications</b>	Western Blot, Simple Western, ELISA, Flow Cytometry, Immunoblotting, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin, Immunoprecipitation, SDS-Page, Proximity Ligation Assay, Chromatin Immunoprecipitation (ChIP), Knockdown Validated, Knockout Validated
<b>Recommended Dilutions</b>	Western Blot, Simple Western, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen, Immunoblotting, Proximity Ligation Assay, SDS-Page, Chromatin Immunoprecipitation (ChIP), Knockout Validated, Knockdown Validated
<b>Application Notes</b>	Optimal dilution of this antibody should be experimentally determined.

## Images

LC3B Antibody [mFluor Violet 450 SE] [NB100-2220MFV450] - Vial of mFluor Violet 450 conjugated antibody. mFluor Violet 450 is optimally excited at 406 nm by the Violet laser (405 nm) and has an emission maximum of 445 nm.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NB100-2220MFV450**

---

NB100-2220PEP	LC3B Antibody Blocking Peptide
NB600-1384PEP	LC3B Antibody Blocking Peptide
H00008878-M01	p62/SQSTM1 Antibody (2C11)
NB100-2331	LC3A Antibody - BSA Free

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-2220MFV450](http://www.novusbio.com/reviews/submit/NB100-2220MFV450)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

