

# Product Datasheet

## WASP [p Ser483, p Ser484] Antibody [Biotin] NB100-2307B

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-2307B](http://www.novusbio.com/NB100-2307B)

Updated 12/12/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-2307B](http://www.novusbio.com/reviews/destination/NB100-2307B)



**NB100-2307B**

WASP [p Ser483, p Ser484] Antibody [Biotin]

| Product Information         |                                                                                                                                                                                                                           |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                    |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                               |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                  |
| Clonality                   | Polyclonal                                                                                                                                                                                                                |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                        |
| Isotype                     | IgG                                                                                                                                                                                                                       |
| Conjugate                   | Biotin                                                                                                                                                                                                                    |
| Purity                      | Immunogen affinity purified                                                                                                                                                                                               |
| Buffer                      | PBS                                                                                                                                                                                                                       |
| Product Description         |                                                                                                                                                                                                                           |
| Host                        | Rabbit                                                                                                                                                                                                                    |
| Gene ID                     | 7454                                                                                                                                                                                                                      |
| Gene Symbol                 | WAS                                                                                                                                                                                                                       |
| Species                     | Human                                                                                                                                                                                                                     |
| Immunogen                   | The immunogen this antibody was made to, is phosphorylated Serine 483, phosphorylated Serine 484 and surrounding residues of human Wiskott-Aldrich syndrome using the numbering given in entry NP_000368.1 (GeneID 7454). |
| Notes                       | [PUBLIC NOTES]                                                                                                                                                                                                            |
| Product Application Details |                                                                                                                                                                                                                           |
| Applications                | Western Blot                                                                                                                                                                                                              |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                              |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                    |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB100-2307B**

---

|                    |                                             |
|--------------------|---------------------------------------------|
| NBP2-24891B        | Rabbit IgG Isotype Control [Biotin]         |
| H00007454-Q01-10ug | Recombinant Human WASP GST (N-Term) Protein |
| 202-IL-010         | IL-2 [Unconjugated]                         |
| NBP2-11140         | WASP Overexpression Lysate                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-2307B](http://www.novusbio.com/reviews/submit/NB100-2307B)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

