

# Product Datasheet

## SMG1 Antibody [DyLight 650]

### NB100-2320C

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-2320C](http://www.novusbio.com/NB100-2320C)

Updated 12/12/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-2320C](http://www.novusbio.com/reviews/destination/NB100-2320C)



**NB100-2320C**

SMG1 Antibody [DyLight 650]

| Product Information         |   |
|-----------------------------|---|
| Unit Size                   | 0.1 ml  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.   |
| Storage                     | Store at 4C in the dark.  |
| Clonality                   | Polyclonal  |
| Preservative                | 0.05% Sodium Azide  |
| Isotype                     | IgG   |
| Conjugate                   | DyLight 650   |
| Purity                      | Immunogen affinity purified   |
| Buffer                      | 50mM Sodium Borate  |
| Product Description         |   |
| Host                        | Rabbit  |
| Gene ID                     | 23049   |
| Gene Symbol                 | SMG1  |
| Species                     | Human   |
| Immunogen                   | The immunogen recognized by this antibody maps to a region between residues 1625 and 1675 of human PI-3-Kinase-related kinase SMG1 using the numbering given in entry NP_055907.3 (GeneID 23049). |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.   |
| Product Application Details |   |
| Applications                | Western Blot, Immunoprecipitation   |
| Recommended Dilutions       | Western Blot, Immunoprecipitation   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.  |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB100-2320C**

---

|                    |   |
|--------------------|---|
| NBP2-24891C        | Rabbit IgG Isotype Control [DyLight 650]    |
| H00023049-Q01-10ug | Recombinant Human SMG1 GST (N-Term) Protein |
| 7754-BH-005/CF     | TGF-beta 1 [Unconjugated]                   |
| NBP1-88913PEP      | SMG1 Recombinant Protein Antigen            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-2320C](http://www.novusbio.com/reviews/submit/NB100-2320C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

