

Product Datasheet

USP4 Antibody - BSA Free NB100-2868

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-2868

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-2868



NB100-2868

USP4 Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)

Product Description


Description	Novus Biologicals Rabbit USP4 Antibody - BSA Free (NB100-2868) is a polyclonal antibody validated for use in IP. Anti-USP4 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	7375
Gene Symbol	USP4
Species	Human
Immunogen	The immunogen recognized by this antibody maps to a region between residue 225 and 275 of human Ubiquitin Specific Protease 4 using the numbering given in entry AAB72237.1 (GeneID 7375).

Product Application Details

Applications	Immunoprecipitation
Recommended Dilutions	Immunoprecipitation 2-5 ug/mg lysate
Application Notes	Not Recommended for Western Blot. Use rabbit anti-USP4 Antibody NB 100-2868.

Images

Immunoprecipitation: USP4 Antibody [NB100-2868] - Detection of Human USP4 by Western Blot of Immunoprecipitates. Samples: Whole cell lysate (1 mg for IP, 20% of IP loaded) from HeLa cells. Antibodies: Affinity purified rabbit anti-USP4 antibody NB100-2868 used for IP at 3 mcg/mg lysate. USP4 was also immunoprecipitated by rabbit anti-USP4 antibody NB100-2869, which recognizes a downstream epitope. Detection: Chemiluminescence with an exposure time of 30 seconds.

**Publications**

Yang Y, Xia F, Hermance N et al. A cytosolic ATM/NEMO/RIP1 complex recruits TAK1 to mediate the NF-kappaB and p38 mitogen-activated protein kinase (MAPK)/MAPK-activated protein 2 responses to DNA damage. Mol Cell Biol 2011-07-01 [PMID: 21606198]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NB100-2868

NBL1-17663	USP4 Overexpression Lysate
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-2868

Earn gift cards/discounts by submitting a publication using this product:

www.novusbio.com/publications

