

Product Datasheet

SH2D1A Antibody NB100-40783

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-40783

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-40783



NB100-40783**SH2D1A Antibody**

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris saline (20 mM Tris pH 7.3, 150 mM NaCl), 0.5% BSA


Product Description	
Host	Goat
Gene ID	4068
Gene Symbol	SH2D1A
Species	Human
Immunogen	Peptide with sequence C-QYPVEKKSSARSTQ corresponding to internal region according to NP_002342.1, NP_001108409.1.

Product Application Details	
Applications	Western Blot, Peptide ELISA
Recommended Dilutions	Western Blot 0.2 - 0.6 ug/ml, Peptide ELISA Detection limit 1:32000
Application Notes	IP: An anonymous customer has reported usage under native conditions. WB: Approx. 16 kDa band observed in human lymph node lysates (calculated MW of 14.2 kDa band according to NP_002342.1). In transfected HEK293 transiently expressing SH2D1A a band of approx. 17 kDa is observed. This band is not observed in the non-transfected HEK293.



Images

Western Blot: SH2D1A Antibody [NB100-40783] - Analysis of SH2D1A using NB100-40783 at 0.2 ug/ml in Human Lymph Node lysate (35ug protein in RIPA buffer).



Western Blot: SH2D1A Antibody [NB100-40783] - HEK293 overexpressing SH2D1A and probed with (mock transfection in first lane).



Publications

Romero X, Benitez D, March S et al. Differential expression of SAP and EAT-2-binding leukocyte cell-surface molecules CD84, CD150 (SLAM), CD229 (Ly9) and CD244 (2B4). Tissue Antigens 2004-08-01 [PMID: 15245368]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-40783

NBL1-15918	SH2D1A Overexpression Lysate
HAF017	Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
HAF109	Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB410-28088-1mg	Goat IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-40783

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



