

# Product Datasheet

## USP24 Antibody [PE] NB100-40830PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-40830PE](http://www.novusbio.com/NB100-40830PE)

Updated 1/10/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-40830PE](http://www.novusbio.com/reviews/destination/NB100-40830PE)



**NB100-40830PE**

USP24 Antibody [PE]

**Product Information**

|                      |   |
|----------------------|---|
| <b>Unit Size</b>     | 0.1 ml  |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C in the dark.  |
| <b>Clonality</b>     | Polyclonal  |
| <b>Preservative</b>  | 0.05% Sodium Azide  |
| <b>Isotype</b>       | IgG   |
| <b>Conjugate</b>     | PE  |
| <b>Purity</b>        | Immunogen affinity purified   |
| <b>Buffer</b>        | PBS   |

**Product Description**

|                    |   |
|--------------------|---|
| <b>Host</b>        | Rabbit  |
| <b>Gene ID</b>     | 23358   |
| <b>Gene Symbol</b> | USP24   |
| <b>Species</b>     | Human   |
| <b>Immunogen</b>   | The immunogen recognized by this antibody maps to a region between residue 2410 and the C-terminus (residue 2460) of human Ubiquitin Specific Peptidase 24 using the numbering given in entry XP_371254.5 (GeneID 23358). |

**Product Application Details**

|                              |  |
|------------------------------|--|
| <b>Applications</b>          | Western Blot, Immunoprecipitation                                      |
| <b>Recommended Dilutions</b> | Western Blot, Immunoprecipitation                                      |
| <b>Application Notes</b>     | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB100-40830PE**

|                   |   |
|-------------------|---|
| NBP2-24983        | Rabbit IgG Isotype Control [PE]                 |
| H00023358-P01-2ug | Recombinant Human USP24 GST (N-Term) Protein    |
| H00023358-T02     | USP24 293T Cell Transient Overexpression Lysate |
| NB600-1160        | UCH-L1/PGP9.5 Antibody (31A3)                   |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-40830PE](http://www.novusbio.com/reviews/submit/NB100-40830PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

