

Product Datasheet

Neurofibromin 1 Antibody [DyLight 350] NB100-418UV

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-418UV

Updated 12/12/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-418UV



NB100-418UV

Neurofibromin 1 Antibody [DyLight 350]

| Product Information | |
|-----------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Polyclonal |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG |
| Conjugate | DyLight 350 |
| Purity | Immunogen affinity purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Rabbit |
| Gene ID | 4763 |
| Gene Symbol | NF1 |
| Species | Human, Mouse, C. elegans |
| Reactivity Notes | C. elegans reactivity reported in scientific literature (PMID: 30929290). |
| Immunogen | The immunogen recognized by this antibody maps to a region between residue 2760 and the C-terminus (residue 2818) of human Neurofibromin 1 using the numbering given in entry NP_000258.1 (GenelD 4763) |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |
| Product Application Details | |
| Applications | Western Blot, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Knockdown Validated |
| Recommended Dilutions | Western Blot, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Knockdown Validated |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-418UV

| | |
|--------------------|--|
| NBP2-24891UV | Rabbit IgG Isotype Control [DyLight 350] |
| NB300-155PEP | Neurofibromin 1 Antibody Blocking Peptide |
| DPSG10 | PSG1 [HRP] |
| H00004763-Q01-10ug | Recombinant Human Neurofibromin 1 GST (N-Term) Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-418UV

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

