

Product Datasheet

NEI3 Antibody - BSA Free NB100-441SS

Unit Size: 0.025 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-441

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-441



NB100-441SS

NEI3 Antibody - BSA Free

Product Information

| | |
|--------------------------------|--|
| Unit Size | 0.025 ml |
| Concentration | 1 mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | 0.02% Sodium Azide |
| Isotype | IgG |
| Purity | Immunogen affinity purified |
| Buffer | PBS |
| Target Molecular Weight | 67.7 kDa |

Product Description


| | |
|--------------------------------|--|
| Host | Rabbit |
| Gene ID | 55247 |
| Gene Symbol | NEIL3 |
| Species | Human |
| Reactivity Notes | Hu |
| Specificity/Sensitivity | This is specific for human NEI3. |
| Immunogen | Residues 164-177 (LRAESEVKKQKGRMLG) of the human NEI3 protein. 85% sequence identity to mouse and rat. |

Product Application Details

| | |
|------------------------------|---|
| Applications | Western Blot |
| Recommended Dilutions | Western Blot 1:500 |
| Application Notes | This antibody may be used for Western blot. No other assays have been tested. Suggested working dilutions: * Western Blot - 1:500 *The investigator should determine the optimal working dilution for a specific application. The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors. |

Images

Western Blot: NEI3 Antibody [NB100-441] - Analysis of NEI3 using IMG-4151 at 1:500 dilution on 50 ug of human testis lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-441SS

| | |
|------------|---|
| NBL1-13581 | NEI3 Overexpression Lysate |
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-441

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

