

# Product Datasheet

## WRN Antibody [mFluor Violet 500 SE] NB100-471MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-471MFV500](http://www.novusbio.com/NB100-471MFV500)

Updated 11/11/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-471MFV500](http://www.novusbio.com/reviews/destination/NB100-471MFV500)



**NB100-471MFV500**

WRN Antibody [mFluor Violet 500 SE]

| Product Information         |                                                                                                                                                                                                                                       |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                           |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                              |
| Clonality                   | Polyclonal                                                                                                                                                                                                                            |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                    |
| Isotype                     | IgG                                                                                                                                                                                                                                   |
| Conjugate                   | mFluor Violet 500 SE                                                                                                                                                                                                                  |
| Purity                      | Immunogen affinity purified                                                                                                                                                                                                           |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                    |
| Product Description         |                                                                                                                                                                                                                                       |
| Host                        | Rabbit                                                                                                                                                                                                                                |
| Gene ID                     | 7486                                                                                                                                                                                                                                  |
| Gene Symbol                 | WRN                                                                                                                                                                                                                                   |
| Species                     | Human                                                                                                                                                                                                                                 |
| Immunogen                   | The immunogen recognized by this antibody maps to a region between residues 400 and 450 of human Werner Syndrome Helicase using the numbering given in SwissProt entry Q14191 (GenelD 7486).                                          |
| Notes                       | mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |
| Product Application Details |                                                                                                                                                                                                                                       |
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP)                                                                                                |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP)                                                                                                |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB100-471MFV500**

---

|                  |                                                   |
|------------------|---------------------------------------------------|
| NBP2-24891MFV500 | Rabbit IgG Isotype Control [mFluor Violet 500 SE] |
| NBP1-87143PEP    | WRN Recombinant Protein Antigen                   |
| NB110-57130      | TRF-2 Antibody - BSA Free                         |
| NB100-309        | ATM Antibody (2C1) - Azide and BSA Free           |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-471MFV500](http://www.novusbio.com/reviews/submit/NB100-471MFV500)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

