

Product Datasheet

NSP 5 alpha 3 alpha Antibody - BSA Free NB100-517

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-517

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-517



NB100-517

NSP 5 alpha 3 alpha Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.1% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)
Target Molecular Weight	80 kDa

Product Description	
Description	Novus Biologicals Rabbit NSP 5 alpha 3 alpha Antibody - BSA Free (NB100-517) is a polyclonal antibody validated for use in WB, ICC/IF and IP. Anti-NSP 5 alpha 3 alpha Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	92521
Gene Symbol	SPECC1
Species	Human
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: mouse (83%), rat (83%).
Immunogen	An internal synthetic peptide made to the human NSP 5 alpha 3 alpha protein sequence (between residues 1-100). [UniProt# Q5M775]


Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:2000, Immunocytochemistry/ Immunofluorescence 1:250, Immunoprecipitation
Application Notes	In Western blot, a band is seen at ~80 kDa. A second band may be seen at ~60 kDa. The specificity of this band is not known.

Images


Immunocytochemistry/Immunofluorescence: NSP 5 alpha 3 alpha Antibody [NB100-517] - Immunostaining of HN30 and MCF-12a cells. MCF-12a NT (DAPI): non-treated MCF-12a DAPI stained, MCF-12a NT (NSP 5a3a-FITC): non-treated MCF-12a FITC stained for NSP 5a3a, MCF-12a pNSP 5a3a (DAPI): HN30 NT (DAPI): non-treated HN30 DAPI stained, HN30 NT (NSP 5a3a-FITC): non-treated HN30 FITC stained for NSP 5a3a, HN30 pNSP 5a3a (DAPI): White arrows indicate apoptotic bodies. Image collected and cropped by CiteAb from the following publication (<https://www.oncotarget.com/article/306/text/>) licensed under a CC-BY license.




Western Blot: NSP 5 alpha 3 alpha Antibody [NB100-517] - Western Blot analysis of total lysates from asynchronous HN30 cells and MCF-12a 3 days post-transfection. NT: non-treated, OB: only buffer, EGFP: only pcDNA3.1/CT-GFP vector, NSP 5a3a/EGFP: pcDNA 3.1/CT-GFP and pcDNA 3.0 NSP 5a3a. Image collected and cropped by CiteAb from the following publication (<https://www.oncotarget.com/article/306/text/>) licensed under a CC-BY license.




Western Blot: NSP 5 alpha 3 alpha Antibody [NB100-517] - Detection of NSP 5 alpha 3 alpha in HeLa whole cell lysate (RIPA) using NB100-517. 0.5ug/ml dilution, 15 second ECL detection.




Immunocytochemistry/Immunofluorescence: NSP 5 alpha 3 alpha Antibody [NB100-517] - The NSP 5 alpha 3 antibody was tested in Hela cells at a 1:250 dilution against Dylight 488 (Green). Alpha tubulin and nuclei were counterstained with Dylight 550 (Red) and DAPI (Blue), respectively.



Immunoprecipitation: NSP 5 alpha 3 alpha Antibody [NB100-517] - Western blot analysis of co-immunoprecipitations of asynchronous HN30 non-treated. NT: non-treated HN30 cells, a. Immunoprecipitation of non-treated HN30 cells using rabbit anti-NSP 5a3a and western blot analysis using rabbit anti-DAXX., c. Immunoprecipitation of non-treated HN30 cells using rabbit anti-NSP 5a3a and western blot analysis using rabbit anti-TRAF2. Image collected and cropped by CiteAb from the following publication (<https://www.oncotarget.com/article/306/text/>) licensed under a CC-BY license.



Western Blot: NSP 5 alpha 3 alpha Antibody [NB100-517] - Western blot analysis of co-immunoprecipitations of asynchronous HN30 non-treated & HN30 cells transfected with pcDNA 3.0 NSP 5a3a & pcDNA 3.1/CT-GFP for 3 days NT: non-treated HN30 cells, TRD: HN30 cells over-expressing NSP 5a3a. a. Immunoprecipitation of non-treated HN30 cells using rabbit anti-NSP 5a3a & western blot analysis using rabbit anti-DAXX., b. Immunoprecipitation of HN30 treated cells using rabbit anti-NSP 5a3a & western blot analysis using rabbit anti-DAXX., c. Immunoprecipitation of non-treated HN30 cells using rabbit anti-NSP 5a3a & western blot analysis using rabbit anti-TRAF2., d. Immunoprecipitation of HN30 treated cells using rabbit anti-NSP 5a3a & western blot analysis using rabbit anti-TRAF2., e. Immunoprecipitation of non-treated HN30 cells using rabbit anti-TRAF2/DAXX & western blot analysis using rabbit anti-NSP 5a3a, & f. Immunoprecipitation of HN30 treated cells using rabbit anti-TRAF2/DAXX & western blot analysis using rabbit anti-NSP 5a3a. Image collected & cropped by CiteAb from the following publication (<https://www.oncotarget.com/lookup/doi/10.18632/oncotarget.306>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

D'Agostino L, Giordano A. A novel dual signaling axis for NSP 5a3a induced apoptosis in head and neck carcinoma *Oncotarget* 2011-12-01 [PMID: 22170762] (IP, WB, ICC/IF, Human)

D'agostino L, Giordano A. NSP 5a3a: A Potential Novel Cancer Target in Head and Neck Carcinoma. *Oncotarget*;1 (6):422-35. 2010-10-01 [PMID: 21311098] (WB, Human)

Procedures

Western Blot protocol for NSP 5 alpha 3 alpha Antibody (NB100-517)

NSP 5 alpha 3 alpha Antibody:

Western Blot Protocol

1. Perform SDS-PAGE (3-8%) on samples to be analyzed, loading 30ug of total protein per lane.
2. Transfer proteins to Nitrocellulose according to the instructions provided by the manufacturer of the transfer apparatus.
3. Stain the blot using ponceau S for 1-2 minutes to access the transfer of proteins onto the nitrocellulose membrane. Rinse the blot in water to remove excess stain and mark the lane locations and locations of molecular weight markers using a pencil.
4. Rinse the blot in TBS for approximately 5 minutes.
5. Block the membrane using 5% non-fat dry milk in TBS for 1 hour.
6. Dilute the rabbit anti-NSP5a3a primary antibody (NB 100-517) in blocking buffer and incubate overnight at 4 degrees Celsius.
7. Wash the membrane in water for 5 minutes and apply the diluted rabbit-IgG HRP-conjugated secondary antibody in blocking buffer (as per manufacturer's instructions) and incubate 1 hour at room temperature.
8. Wash the blot in TBS containing 0.05-0.1% Tween-20 for 10-20 minutes.
9. Wash the blot in type I water for an additional 10-20 minutes (this step can be repeated as required to reduce background).
10. Apply the detection reagent of choice in accordance with the manufacturer's instructions (Amersham's ECL is the standard reagent used at Novus Biologicals).

Note: Tween-20 can be added to the blocking buffer at a final concentration of 0.05-0.2%, provided it does not interfere with antibody-antigen binding.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-517

NB800-PC1	HeLa Whole Cell Lysate
NB100-517PEP	NSP 5 alpha 3 alpha Antibody Blocking Peptide
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-517

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



