

Product Datasheet

TRAPPC9 Antibody NB100-55735

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-55735

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-55735



NB100-55735**TRAPPC9 Antibody**

Product Information	
Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Sodium Azide
Isotype	IgG
Purity	Unpurified
Buffer	Whole Rabbit antisera with 50% Glycerol

Product Description	
Host	Rabbit
Gene ID	83696
Gene Symbol	TRAPPC9
Species	Human, Mouse, Bovine, Canine, Chicken, Primate
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: 100% homologous in human (isoform CRA_A and CRA_b), chimpanzee, monkey, mouse, opossum, rat and dog; xenopus (94%).
Immunogen	Amino acids 282-299 (KKDFVGLDTSRHYKKRC) of human NIBP was used as the immunogen.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:100-1:2000. Use reported in scientific literature (Zahoor et al (2010)), Immunocytochemistry/ Immunofluorescence reported in scientific literature (Zahoor et al (2010)), Immunoprecipitation 1:10-1:500. Use reported in scientific literature (Zahoor et al (2010))
Application Notes	Confocal immunofluorescence: Cells were fixed with 4% paraformaldehyde in PBS for 10 min and then permeabilized with 0.5 % Triton X-100 for 20 min at room temperature (RT). The cells were blocked with PBS containing 10 % BSA and probed with NIBP pAB for 1 hour at RT. After washing with PBS, cells were incubated with Alexa conjugated goat anti-mouse antibodies. See Zahoor et al for details (2010).

Publications

Zahoor MA, Yamane D, Mohamed YM et al. Bovine viral diarrhea virus non-structural protein 5A interacts with NIK- and IKKbeta-binding protein. J Gen Virol. 2010-08-01 [PMID: 20444997] (WB)

Details:

IF (Fig 4): NIBP-FLAG transiently transfected LB9.K (bovine epithelial) cell line. Note: NIBP pAB was transfected validated in LB9.K cells by immunofluorescence. IP & IP/WB (Fig 2): Endogenous NIBP in LB9.K cells infected with bovine viruses was immunopre





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-55735

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00083696-P01-10ug	Recombinant Human TRAPPC9 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-55735

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

