

Product Datasheet

CXCR3 Antibody - BSA Free

NB100-56404

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 11

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-56404

Updated 2/21/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-56404



NB100-56404


CXCR3 Antibody - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	41 kDa
Product Description	
Host	Rabbit
Gene ID	2833
Gene Symbol	CXCR3
Species	Human, Mouse
Immunogen	This antibody was developed against a synthetic peptide corresponding to amino acids 290-303 of human CXCR3.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin, Proximity Ligation Assay
Recommended Dilutions	Western Blot 1 - 3 ug/mL, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence reported in scientific literature (PMID 31887215), Immunohistochemistry-Paraffin 1:10 - 1:500, Immunohistochemistry-Frozen reported in scientific literature (PMID 19234157), Proximity Ligation Assay reported in scientific literature (PMID 34433033)
Application Notes	The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.




Images


Western Blot: CXCR3 Antibody [NB100-56404] - HEK293T cells.
Western blot image submitted by a verified customer review.




Immunocytochemistry/Immunofluorescence: CXCR3 Antibody [NB100-56404] - Brain microvascular endothelial cells isolated from WT and CXCR3 KO mice were stained with anti-CXCR3 (green) and anti-ZO-1 (blue) antibodies to analyze CXCR3 expression by immunofluorescence. Image from verified customer review.



Immunohistochemistry-Paraffin: CXCR3 Antibody [NB100-56404] - IHC analysis of a FFPE human spleen using 1:25 dilution of CXCR3 antibody on a Bond Rx autostainer (Leica Biosystems). The assay involved 20 minutes of HIER using 10mM sodium citrate buffer (pH 6.0) and endogenous peroxidase quenching with peroxide block. The sections were incubated with primary antibody for 30 minutes and Bond Polymer Refine Detection (Leica Biosystems) with DAB was used for signal development followed by counterstaining with hematoxylin. Whole slide scanning and capturing of representative images (40X) was performed using Aperio AT2 (Leica Biosystems). Staining was performed by Histowiz.



Western Blot: CXCR3 Antibody [NB100-56404] - Analysis of CXCR3 in human pancreas lysate using CXCR3 antibody at 2 ug/mL. Goat anti-rabbit Ig HRP secondary antibody and PicoTect ECL substrate solution were used for this test.



Publications

Mendes N, Zanesco A, Aguiar C et al. CXCR3-expressing myeloid cells recruited to the hypothalamus protect against diet-induced body mass gain and metabolic dysfunction *Elife* 2024-11-13 [PMID: 39535032]

Lin Y, Zhou HC, Chen N et al. Unveiling the improved targeting migration of mesenchymal stem cells with CXCR3 chemokine receptor 3-modification using intravital NIR-II photoacoustic imaging *Journal of Nanobiotechnology* 2022-12-01 [PMID: 35764961]

Alluri SR, Higashi Y, Berendzen A et al. Synthesis and preclinical evaluation of a novel fluorine-18 labeled small-molecule PET radiotracer for imaging of CXCR3 receptor in mouse models of atherosclerosis *EJNMMI Research* 2023-07-13 [PMID: 37438543] (Immunohistochemistry)

Alluri SR, Higashi Y, Berendzen A et al. Synthesis and preclinical evaluation of a novel fluorine-18 labeled small-molecule PET radiotracer for imaging of CXCR3 receptor in mouse models of atherosclerosis *Research square* 2023-02-23 [PMID: 36865232] (IHC, Mouse)

Zhu L, Gao H, Huang S et al. Combination of STING agonist and CXCR3 antagonist disrupts immune tolerance to overcome anti-PD-L1 resistance in lung adenocarcinoma under oxidative stress *Gene* 2022-10-19 [PMID: 36272651] (IF/IHC, Mouse)

Park JY, Yoo KD, Bae E et al. Blockade of STAT3 signaling alleviates progression of acute kidney injury-to-chronic kidney disease through anti-apoptosis *American journal of physiology. Renal physiology* 2022-03-21 [PMID: 35311382] (IF/IHC)

Javeed N, Her TK, Brown MR Et al. Pro-inflammatory beta cell small extracellular vesicles induce beta cell failure through activation of the CXCL10/CXCR3 axis in diabetes *Cell reports* 2021-08-24 [PMID: 34433033] (PLA, Mouse)

Petrisko TJ, Bloemer J, Pinky PD et al. Neuronal CXCL10/CXCR3 Axis Mediates the Induction of Cerebral Hyperexcitability by Peripheral Viral Challenge *Front Neurosci* 2020-03-24 [PMID: 32265633] (ICC/IF, Mouse)

Alvarez-Carbonell D, Ye F, Ramanath N et al. Cross-talk between microglia and neurons regulates HIV latency *PLoS Pathog.* 2019-12-01 [PMID: 31887215] (ICC/IF, Human)

Suttorp CM, Cremers NA, van Rheden R et al. Chemokine Signaling during Midline Epithelial Seam Disintegration Facilitates Palatal Fusion *Front Cell Dev Biol.* 2017-10-30 [PMID: 29164113] (IHC-P, Mouse)

Details:

Used CXCR3 IHC antibody from Novus to quantify CXCR3 in palatal shelve mesenchyme

Brambilla R, Persaud T, Hu XC et al. Transgenic Inhibition of Astroglial NF- κ B Improves Functional Outcome in Experimental Autoimmune Encephalomyelitis by Suppressing Chronic Central Nervous System Inflammation . *J Immunol* 182 2628-2640. 2009-01-01 [PMID: 19234157]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-56404

NBL1-09631	CXCR3 Overexpression Lysate
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-56404

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

