

Product Datasheet

TLR4 Antibody (76B357.1) [Janelia Fluor® 646] NB100-56566JF646

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-56566JF646

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-56566JF646



NB100-56566JF646

TLR4 Antibody (76B357.1) [Janelia Fluor® 646]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	76B357.1
Preservative	0.05% Sodium Azide
Isotype	IgG2b Kappa
Conjugate	Janelia Fluor 646
Purity	Protein G purified
Buffer	50mM Sodium Borate
Target Molecular Weight	95.7 kDa
Product Description	
Host	Mouse
Gene ID	7099
Gene Symbol	TLR4
Species	Human, Mouse, Rat, Porcine, Bovine, Mammal
Immunogen	This TLR4 Antibody (76B357.1) was developed against a portion of amino acids 100-200 of human TLR4 (NP_612564).
Notes	Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.
Product Application Details	
Applications	Western Blot, Dot Blot, ELISA, Flow Cytometry, Flow (Cell Surface), Flow (Intracellular), Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin, In vitro assay, SDS-Page, Block/Neutralize, Chromatin Immunoprecipitation (ChIP), CyTOF-ready, Dual RNAscope ISH-IHC, ELISA Capture (Matched Antibody Pair), Knockdown Validated, Knockout Validated
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen, In vitro assay, Dot Blot, SDS-Page, Flow (Cell Surface), Flow (Intracellular), ELISA Capture (Matched Antibody Pair), Chromatin Immunoprecipitation (ChIP), CyTOF-ready, Knockout Validated, Knockdown Validated, Block/Neutralize, Dual RNAscope ISH-IHC
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-56566JF646

NBP1-43317JF646	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b) [Janelia Fluor 646]
NBP2-24821PEP	TLR4 Antibody Blocking Peptide
210-TA-005	TNF-alpha [Unconjugated]
NBP2-26244	TLR4 Inhibitor Peptide Set

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-56566JF646

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

