

Product Datasheet

NCOR2 Antibody NB100-58826

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-58826

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-58826



NB100-58826

NCOR2 Antibody


Product Information	
Unit Size	0.1 ml
Concentration	0.2 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	TBS and 0.1% BSA

Product Description	
Description	This antibody can be used as the primary antibody in a PLA assay with the following as complementing antibodies:NB500-160, NB100-381, NB200-347, NBP1-06527
Host	Rabbit
Gene ID	9612
Gene Symbol	NCOR2
Species	Human
Immunogen	The immunogen recognized by this antibody maps to a region between residue 1525 and 1575 of human silencing mediator of retinoic acid and thyroid hormone receptor alpha using the numbering given in entry AAD20946.1 (GeneID 9612).
Notes	This antibody can be used as the primary antibody in a PLA assay with the following as complementing antibodies:NB500-160, NB100-381, NB200-347, NBP1-06527

Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:2000-1:10000, Immunoprecipitation 2-5 ug/mg lysate

Images

Western Blot: NCOR2 Antibody [NB100-58826] - Detection of Human SMRT on HeLa whole cell lysate . SMRT was also immunoprecipitated by rabbit anti-SMRT antibody NB100-58827.



Publications

Abe Y, Kofman ER, Almeida M et al. RANK ligand converts the NCoR/HDAC3 co-repressor to a PGC1?- and RNA-dependent co-activator of osteoclast gene expression *Molecular cell* 2023-09-14 [PMID: 37751740]

Zhou R, Li X, Hu G et al. miR-16 targets transcriptional corepressor SMRT and modulates NF-kappaB-regulated transactivation of interleukin-8 gene. *PLoS One* 2012-01-01 [PMID: 22292036]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-58826

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00009612-Q01-10ug	Recombinant Human NCOR2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-58826

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

