

Product Datasheet

TRIM37 Antibody - BSA Free NB100-59820

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-59820

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-59820



NB100-59820

TRIM37 Antibody - BSA Free

Product Information

| | |
|----------------------|---------------------------------------|
| Unit Size | 100 ul |
| Concentration | 1.0 mg/ml |
| Storage | Store at 4C. Do not freeze. |
| Clonality | Polyclonal |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG |
| Purity | Immunogen affinity purified |
| Buffer | Tris-Citrate/Phosphate (pH 7.0 - 8.0) |

Product Description


| | |
|--------------------|--|
| Host | Rabbit |
| Gene ID | 4591 |
| Gene Symbol | TRIM37 |
| Species | Human |
| Immunogen | The immunogen recognized by this antibody maps to a region between residue 850 and 900 of human tripartite motif-containing 37 using the numbering given in entry NP_056109.1 (GenelD 4591). |

Product Application Details

| | |
|------------------------------|---|
| Applications | Western Blot, Immunoprecipitation |
| Recommended Dilutions | Western Blot, Immunoprecipitation 2- 5 ug/mg lysate |
| Application Notes | Use in western blot reported in scientific literature (PMID 24317724) |

Images

Immunoprecipitation: TRIM37 Antibody [NB100-59820] - Detection of Human TRIM37 on HeLa whole cell lysate using NB100-59820. TRIM37 was also immunoprecipitated by rabbit anti-TRIM37 antibody NB100-59821. For blotting immunoprecipitated TRIM37 NB100-59821 was used at 1.0 ug/ml.

**Publications**

Tabah AA, Tardif K, Mansky LM. Anti-HIV-1 Activity of Trim 37. J Gen Virol 2013-12-06 [PMID: 24317724] (WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-59820

| | |
|------------|---|
| NBL1-17293 | TRIM37 Overexpression Lysate |
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-59820

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

