

Product Datasheet

KAT3B/p300 Antibody (RW105) [FITC] NB100-616F-0.1ml

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-616F

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-616F



NB100-616F-0.1ml

KAT3B/p300 Antibody (RW105) [FITC]

| Product Information | |
|---------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | RW105 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 Kappa |
| Conjugate | FITC |
| Purity | Protein G purified |
| Buffer | PBS |

| Product Description | |
|-------------------------|--|
| Host | Mouse |
| Gene ID | 2033 |
| Gene Symbol | EP300 |
| Species | Human, Mouse, Rat, Mustelid, Primate |
| Reactivity Notes | Mink. |
| Specificity/Sensitivity | This is specific for p300 protein. This recognizes residues 1921-2023. |
| Immunogen | Fusion protein containing residues 1572-2371 of human KAT3B/p300. [UniProt# Q0947] |

| Product Application Details | |
|-----------------------------|--|
| Applications | Western Blot, Electron Microscopy, Flow Cytometry, Flow (Intracellular), Immunocytochemistry/ Immunofluorescence, In vitro assay, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP), Knockdown Validated |
| Recommended Dilutions | Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, In vitro assay, Electron Microscopy, Flow (Intracellular), Chromatin Immunoprecipitation (ChIP), Knockdown Validated |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-616F-0.1ml

| | |
|--------------------|--|
| NBP1-43319F | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1) [FITC] |
| NB100-616PE-0.1ml | KAT3B/p300 Antibody (RW105) [PE] |
| H00002033-Q01-10ug | Recombinant Human KAT3B/p300 GST (N-Term) Protein |
| NBP2-29661-1Kit | Human, Mouse, Rat RelA/NFkB p65 ELISA Kit (Colorimetric) |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-616F

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

